



Reusable Silicone Lab Mat



Product Highlights

- Made from durable silicone
- Versatile! Use on the bench, in cold rooms, sterile suites and areas where high sterility is a requirement
- Anti-skid surface
- Noise dampening! Reduces vibration sound from common worktop equipment
- Protects your bench from hot objects up to 200°C
- Chemically resistant and chemically inert
- Can be rolled up for easy storage
- Autoclavable

Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



General Data

Art. No.	E3860-0061
Pack Size	1 Piece (1 Pack × 1 Piece)
Dimensions (W × D × H)	60 cm × 35 cm × 2 cm
Dimensions (W × D)	60 cm × 35 cm
External width	60 cm
Autoclavable	Yes
Material	Silicone
Color	One side solid yellow, other side blue / grey design
Temperature use	-50 °C to 200 °C
Dishwasher safe	Yes
Reusable	Yes
Reversible	Yes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



More informations about Reusable Silicone Lab Mat

Versatile reusable mat

A reusable mat for your lab bench to help keep it clean and safe from stains and spills. Made from durable silicone, the lab mat creates a stain-resistant washable working surface. After use, the mat can easily be cleaned using any standard disinfectant or simply by wiping if no contaminants are present. It can also be put into the dishwasher and it can be autoclaved.

The anti-skid mat offers a large and safe working area when using harsh chemicals and delicate glassware. The lipped edge of the mat helps contains spills and prevents liquids from running onto the floor.