




## 96-Well Aluminium Plate



### Product Highlights

- 96-well aluminium plate
- For use with 0.2 ml PCR tubes or one 96-well plate
- Use on the bench to aid sample preparation or as a cooled storage block
- Cool the block down in the freezer before using as a chiller rack on the bench

## Product Applications

Sample storage

PCR sample preparation

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## General Data

Art. No.	N0000-0096
Pack Size	1 Piece (1 Rack × 1 Piece)
Dimensions (W × D × H)	114.3 mm × 76.2 mm × 21.6 mm
Dimensions (W × D)	114.3 mm × 76.2 mm
Rack Type	Chiller Rack
Material	Aluminium
Color	Silver
Places	96

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)