



Filtration Funnels (sterile)



Product Highlights

- Pore sizes 0.22 and 0.45 µm
- High speed filtration
- Compatible with all common bottles (GL 45 compatible)
- Temperature range: -20 °C to +50 °C
- Sterile SAL10⁻⁶ (ISO1137) / 3-year shelf life
- Certified RNase, DNase free and non-pyrogenic
- Hose connector designed to fit multiple hose diameters

Product Applications

Sterile filtration
Media preparation

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



General Data

Art. No.	See variations
Pack Size	24 Pcs. (1 Box × 24 Pcs.)
Sterile	Yes
Autoclavable	No
DNase free	Yes
Free of endotoxins	Yes
Free of heavy metals	Yes
Pyrogen free	Yes
RNase free	Yes
Non-cytotoxic	Yes
System	Funnel Only

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



More informations about Filtration Funnels (sterile)













Bottle Top Filtration Funnels with PES mebrane for fast filtration. The clever hose connector is designed to fit multiple hose diameters. Each product is manufactured in a Class 100,000 cleanroom environment and is printed with the batch number and expiry date for easy traceability.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk







All Variations

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	Bottle Top Filtration Unit, Funnel Only, 0.22 µm, 150 ml (sterile) Pore size: 0.22 µm Volume: 150 ml	24 Pcs. (1 Box × 24 Pcs.)	CC8222-3226
	 Create PDF Data Sheet		
	Bottle Top Filtration Unit, Funnel Only, 0.22 µm, 250 ml (sterile) Pore size: 0.22 µm Volume: 250 ml	24 Pcs. (1 Box × 24 Pcs.)	CC8222-4226
	 Create PDF Data Sheet		
	Bottle Top Filtration Unit, Funnel Only, 0.22 µm, 500 ml (sterile) Pore size: 0.22 µm Volume: 500 ml	24 Pcs. (1 Box × 24 Pcs.)	CC8222-5226
	 Create PDF Data Sheet		
	Bottle Top Filtration Unit, Funnel Only, 0.22 µm, 1000 ml (sterile) Pore size: 0.22 µm Volume: 1000 ml	24 Pcs. (1 Box × 24 Pcs.)	CC8222-6226
	 Create PDF Data Sheet		
	Bottle Top Filtration Unit, Funnel Only, 0.45 µm, 150 ml (sterile) Pore size: 0.45 µm Volume: 150 ml	24 Pcs. (1 Box × 24 Pcs.)	CC8222-3456
	 Create PDF Data Sheet		
	Bottle Top Filtration Unit, Funnel Only, 0.45 µm, 250 ml (sterile) Pore size: 0.45 µm Volume: 250 ml	24 Pcs. (1 Box × 24 Pcs.)	CC8222-4456
	 Create PDF Data Sheet		

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



PRODUCT NAME		PACKAGING SIZE	ART. NO.
	Bottle Top Filtration Unit, Funnel Only, 0.45 µm, 500 ml (sterile) Pore size: 0.45 µm Volume: 500 ml	24 Pcs. (1 Box × 24 Pcs.)	CC8222-5456
	 Create PDF Data Sheet		
	Bottle Top Filtration Unit, Funnel Only, 0.45 µm, 1000 ml (sterile) Pore size: 0.45 µm Volume: 1000 ml	24 Pcs. (1 Box × 24 Pcs.)	CC8222-6456
	 Create PDF Data Sheet		

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk