



Media Bottles (sterile)



Product Highlights

- Volume 150, 250, 500 and 1000 ml
- Compatible with all common Bottle Top Filters (GL 45 compatible)
- Wide base for stability
- Temperature range: -20 °C to +50 °C
- Sterile SAL10⁻⁶ (ISO1137) / 3-year shelf life
- Certified RNase, DNase free and non-pyrogenic
- Wide bottle top for fast and efficient use
- Ergonomically designed for easier handling to help prevent spills.

Product Applications

Sterile filtration Media preparation Media storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD

5 Tanners Drive MK14 5BU Milton Keynes Email: infoline@starlab.co.uk







General Data

Art. No.	See variations		
Pack Size	24 Pcs. (1 Bag × 24 Pcs.)		
Sterile	Yes		
Autoclavable	No		
Material	Polystyrene		
DNase free	Yes		
Free of endotoxins	Yes		
Free of heavy metals	Yes		
Pyrogen free	Yes		
RNase free	Yes		
Non-cytotoxic	Yes		
System	Media Bottle		

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD

5 Tanners Drive MK14 5BU Milton Keynes Email: infoline@starlab.co.uk







More informations about Media Bottles (sterile)

OPEN WIDE!

Media Bottles with a wide base for stability and a wide, easily accessible bottle top for efficient use. The reservoir bottle caps have large knurls for easy loosening or tightening. Each product is manufactured in a Class 100,000 cleanroom environment and is printed with the batch number and expiry date for easy traceability.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD 5 Tanners Drive

MK14 5BU Milton Keynes Email: infoline@starlab.co.uk







All Variations

PRODUCT NAM	E	PACKAGING SIZE	ART. NO.
	Media Bottle, PS, 150 ml (sterile) Volume: 150 ml Create PDF Data Sheet	24 Pcs. (1 Box × 24 Pcs.)	CC8221-3000
ā	Media Bottle, PS, 250 ml (sterile) Volume: 250 ml Create PDF Data Sheet	24 Pcs. (1 Box × 24 Pcs.)	CC8221-4000
ā	Media Bottle, PS, 500 ml (sterile) Volume: 500 ml Create PDF Data Sheet	24 Pcs. (1 Box × 24 Pcs.)	CC8221-5000
Û	Media Bottle, PS, 1000 ml (sterile) Volume: 1000 ml Create PDF Data Sheet	24 Pcs. (1 Box × 24 Pcs.)	CC8221-6000

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.



5 Tanners Drive MK14 5BU Milton Keynes Email: infoline@starlab.co.uk

