




## Polyolefin Sealing Film for qPCR, Self-Adhesive



### Product Highlights

- Suitable for Real-Time / qPCR
- High optical clarity
- Low autofluorescence
- Low tack-to-touch adhesive
- DMSO, acetonitrile and methanol resistant
- Ideal for dry samples (eg. soil, seeds, insects)
- Temperature range: -80 °C to +105 °C
- Perforated end tabs for easy application and removal
- Certified RNase & DNase free

## Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data

Art. No.	E2796-9895
Overall Dimensions	76.2 mm × 133.4 mm
Working Dimensions	76.2 mm × 133.4 mm
Sterile	No
Pack Size	100 Pcs. (1 Box × 100 Pcs.)
Material	Polyolefin
DNase free	Yes
Free of heavy metals	No
Pyrogen free	No
RNase free	Yes
Color	Clear
Plate type compatibility	Polypropylene, Polystyrene, Polycarbonate
Max. temp.	105 °C
Temperature range	-80 to 105 °C
Temperature min. (SL)	-80 °C
Breathable	no
Easy to Pierce with Needle or Metal Probe	no
Easy to Pierce with Pipette Tip	no
Low 'Tack to Touch' Adhesive	yes
Low Auto-Fluorescence	yes
Opaque	no
Optically clear	yes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



Resistant to DMSO	yes
Seal	Self-adhesive
Short-Term Storage and Incubation	yes
Suitable for PCR	yes
Suitable for Real-Time PCR (qPCR)	yes
Superior Optical Clarity	yes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## More informations about Polyolefin Sealing Film for qPCR, Self-Adhesive

OPTICALLY CLEAR SEAL WITH PRESSURE-ACTIVATED ADHESIVE

Self-adhesive plate seal with superior optical clarity and low fluorescence. The low tack-to-touch adhesive is not tacky until pressure is applied to the sealing surface, which then creates a secure seal around both the wells and the edges of the plate. The seal features end tabs with perforations for easy application and removal. Suitable for qPCR, and for dry samples such as soil, seeds or insects. Use for qPCR, long-term storage, protein crystallography and HPLC.



A roller or seal applicator is recommended for use with this seal to ensure consistent and adequate pressure.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## Accessories

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">Microtitre Sealing Brayer</a>	1 Piece (1 Pack × 1 Piece)	E9127-2940
 <a href="#">Plate Seal Applicator</a>	5 Pcs. (1 Pack × 5 Pcs.)	I2928-7355

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)