



## Starlab Minicentrifuge



### Product Highlights

- 8-place rotor for 1.5 or 2.0 ml tubes, including screw cap tubes
- 32-place rotor for 0.2 ml PCR strips or tubes
- Magnetic lid lock for easy, one-handed operation
- Quick start! Spinning starts on lid closure
- Smooth, quiet operation
- One-touch stop for fast braking
- Automatic lid opening on stop
- Quick and easy, tool-free rotor exchange
- 2-year warranty

## Product Applications

Sample preparation

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data

Art. No.	S8030-5000
Dimensions (W × D × H)	160 mm × 130 mm × 112 mm
Dimensions (W × D)	160 mm × 130 mm
Weight	1.2 kg
Pack Size	1 Piece (1 Box × 1 Piece)
Power supply	100 – 240 V, 50/60 Hz
Operating temperature	2 °C – 40 °C
RCF	2,500 × g

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## More informations about Starlab Minicentrifuge

### **SMOOTH OPERATOR**

With a sleek, curved design that will look fantastic on your bench, the Starlab Minicentrifuge is a compact and solid unit that delivers a quiet, smooth operation.

### **AHEAD OF THE CURVE**

The Starlab Minicentrifuge has a compact footprint that makes it the perfect, personal bench partner for your quick spin downs. It's so easy to use! Just load your samples, close the lid and you're spinning! To end the cycle, touch the Stop button for super fast braking and release of the lid lock. We even thought about your samples. Should the power be cut during spinning, the lid will automatically unlock so your precious samples are not trapped inside. We really have thought of everything!

The Minicentrifuge is supplied with a Power Adapter and four exchangeable plugs, 1 x Tube Rotor, 1 x PCR Strip Rotor, and 8 x 0.5 ml and 8 x 0.2 ml Adapters for use with the Tube Rotor.



## Accessories

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">Starlab Minicentrifuge Tube Rotor Adapters for 0.2 ml tubes</a>	8 Pcs. (1 Box × 8 Pcs.)	S8030-5014
 <a href="#">Starlab Minicentrifuge Tube Rotor Adapters for 0.5 ml tubes</a>	8 Pcs. (1 Box × 8 Pcs.)	S8030-5013
 <a href="#">Power Adapter for Starlab Minicentrifuge or Vortex Mixer</a>	1 Piece (1 Box × 1 Piece)	S8030-5012
 <a href="#">PCR Strip Rotor for Starlab Minicentrifuge</a>	1 Piece (1 Box × 1 Piece)	S8030-5011
 <a href="#">Tube rotor for Starlab Minicentrifuge</a>	1 Piece (1 Box × 1 Piece)	S8030-5010

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)