



## Starlab Minicentrifuge Tube Rotor Adapters for 0.5 ml tubes



### Product Highlights

- Converts the Tube Rotor for use with 0.5 ml tubes
- Pack of 8 adapters
- For use with Starlab Minicentrifuge Tube Rotor only

## Product Applications

Sample

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data

|                     |                         |
|---------------------|-------------------------|
| Art. No.            | S8030-5013              |
| Alternative Product | S8030-5014              |
| Pack Size           | 8 Pcs. (1 Box × 8 Pcs.) |
| RCF                 | 2,500 × g               |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)