



## TubeOne® Microcentrifuge Tubes, 2.0 ml



### Product Highlights

- Conservatively rated to 20,000 x g, but will perform to higher specification under optimal conditions
- Tubes stay sealed during boiling, freezing or centrifuging
- Certified RNase, DNase, DNA and Pyrogen free
- Caps feature extended sealing zone for secure sealing
- Easy to read graduations
- Frosted writing area on side
- Autoclavable
- Manufactured from 100 % pure virgin polypropylene
- Supplied in resealable bags

## Product Applications

Sample preparation

Centrifugation

Sample boiling or freezing

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## General Data

|                       |                |
|-----------------------|----------------|
| Art. No.              | See variations |
| Sterile               | No             |
| Autoclavable          | Yes            |
| Volume                | 2.0 ml         |
| Material              | Polypropylene  |
| DNA free              | Yes            |
| DNase free            | Yes            |
| Free of endotoxins    | Yes            |
| Free of heavy metals  | Yes            |
| Pyrogen free          | Yes            |
| RNase free            | Yes            |
| Color                 | Various        |
| Operating temperature | -80°C to 121°C |
| g-Force max.          | 20,000 x g     |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)











## More informations about TubeOne® Microcentrifuge Tubes, 2.0 ml

Starlab TubeOne® 2.0 ml microcentrifuge tubes protect your samples from evaporation, leakage, accidental opening or breakage during boiling, freezing centrifugation and other critical protocol steps. Made from 100% pure polypropylene, TubeOne® tubes contain no processing additives such as slip agents, biocides or antistatic agents - your samples remain protected for potential bioactive contaminants.



## All Variations

| PRODUCT NAME   | PACKAGING SIZE                     | ART. NO.     |
|--|------------------------------------|--------------|
| <br><a href="#">2.0ml TubeOne® Microcentrifuge Tubes, Natural</a><br><b>Color:</b> Natural<br><b>Pack Size:</b> 500 Tubes (1 Bag × 500 Tubes)<br> <a href="#">Create PDF Data Sheet</a>                          | 500 Tubes (1 Bag × 500 Tubes)      | S1620-2700   |
| <br><a href="#">2.0ml TubeOne® Microcentrifuge Tubes, Natural</a><br><b>Color:</b> Natural<br><b>Pack Size:</b> 15,000 Tubes (30 Bags × 500 Tubes)<br> <a href="#">Create PDF Data Sheet</a>                     | 15,000 Tubes (30 Bags × 500 Tubes) | S1620-2700-C |
| <br><a href="#">2.0ml TubeOne® Microcentrifuge Tubes, Mixed</a><br><b>Color:</b> Mixed<br><b>Pack Size:</b> 500 Tubes (5 Bags × 100 Tubes)<br> <a href="#">Create PDF Data Sheet</a>                          | 500 Tubes (5 Bags × 100 Tubes)     | S1620-2799   |
| <br><a href="#">2.0ml TubeOne® Microcentrifuge Tubes, Mixed, Case (non-sterile)</a><br><b>Color:</b> Mixed<br><b>Pack Size:</b> 12,000 Tubes (24 Bags × 500 Tubes)<br> <a href="#">Create PDF Data Sheet</a> | 12,000 Tubes (24 Bags × 500 Tubes) | S1620-2799-C |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## Accessories

| PRODUCT NAME  |  | PACKAGING SIZE              | ART. NO.   |
|---|--|-----------------------------|------------|
|  | <a href="#">Cap-Locks for 1.5 / 2.0 ml Microcentrifuge Tubes</a> | 100 Pcs. (1 Bag × 100 Pcs.) | I1415-1508 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)