




## Storage box 5.0ml tubes, blue



### Product Highlights

- 25-place storage box for 5.0 ml Centrifuge Tubes
- Transparent lid with printed reference grid for easier sample identification
- One-way keyed lid prevents sample mix-up
- Air vents minimize condensation
- Base with drain holes
- Stackable
- Temperature use: -90°C to 80°C
- Dimensions (WxDxH): 133 x 133 x 67 mm (with lid)
- Cover made of polystyrene / Base made of high impact polystyrene

## Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data

Art. No.	E2325-5021
Dimensions (W × D × H)	133 mm × 133 mm × 67 mm
Dimensions (W × D)	133 mm × 133 mm
Autoclavable	No
Pack Size	3 Pcs. (1 Pack × 3 Pcs.)
Material	Polystyrene
Color	Blue

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## More informations about Storage box 5.0ml tubes, blue

25-place stackable, storage box for your 5.0 ml Centrifuge Tubes. The transparent lid features a printed reference grid for easy sample identification Air vents minimise condensation and drain holes in the base aides lid removal following freezing. Suitable for use down to  $-90^{\circ}\text{C}$ .

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)