



Centrifuge tube, 5ml, screw tube only, conical, natural



Product Highlights

- RNase, DNase, DNA & Pyrogen free
- Temperature use: -196°C to 121°C
- Autoclavable
- Centrifugation rate: 25,000 x g
- Molded graduations at 0.5 ml increments
- Etched writing surface on side wall
- Dimensions (HxØ): 62.5 x 15.8 mm
- Made of ultra-clear polypropylene
- Caps available separately

Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



General Data

Art. No.	E1450-5040
Pack Size	200 Tubes (2 Bags × 100 Tubes)
Sterile	No
Volume	5.0 ml
Cap style	Screw Cap
Material	Polypropylene
DNA free	Yes
DNase free	Yes
Pyrogen free	Yes
RNase free	Yes
Color	Natural
g-Force max.	25.000 x g

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



More informations about Centrifuge tube, 5ml, screw tube only, conical, natural











Do you need to prepare samples up to 5.0 ml but don't want to use 15 ml tubes? Save on plastic waste, with our 5.0 ml Centrifuge Tubes!

- Ideal for biological or chemical solution preparation
- Easier sample access lowers the risk of contamination
- Saves on space and plastic, compared to 15 ml tubes
- Suitable for microbial incubation and pellet preparation
- Choice of cap styles available - flip top or screw cap
- Sterile options available
- Amber tube for light sensitive reactions
- Various accesories available

No processing additives such as slip agents, biocides, plasticizers or mold release agents are used during production.



All Variations





PRODUCT NAME		PACKAGING SIZE	ART. NO.
	Centrifuge tube, 5ml, flip cap, conical, natural Cap style: Flip Cap Sterile: No Color: Natural  Create PDF Data Sheet	200 Tubes (2 Bags × 100 Tubes)	E1450-5000
	Centrifuge tube, 5ml, flip cap, conical, amber Cap style: Flip Cap Sterile: No Color: Amber  Create PDF Data Sheet	200 Tubes (2 Bags × 100 Tubes)	E1450-5007
	Centrifuge tube, 5ml, flip cap, conical, natural, sterile Cap style: Flip Cap Sterile: Yes Color: Natural  Create PDF Data Sheet	200 Tubes (10 Bags × 20 Tubes)	E1450-5010
	Centrifuge tube, 5ml, screw tube only, conical, natural Cap style: Screw Cap Sterile: No Color: Natural  Create PDF Data Sheet	200 Tubes (2 Bags × 100 Tubes)	E1450-5040
	Centrifuge tube, 5.0ml, with standard screw cap, conical, natural white, sterile Cap style: Screw Cap Sterile: Yes Color: Natural  Create PDF Data Sheet	200 Tubes (10 Bags × 20 Tubes)	E1450-5050

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



Accessories

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	Rack 5.0ml tubes, blue	5 Pcs. (1 Pack × 5 Pcs.)	E2325-5011
	Universal centrifuge adapter, for 5.0 ml tubes, for 15 ml rotor bores	1 Pack × 8 Pcs.	E1450-9011
	Universal centrifuge adapter, for 5.0 ml tubes, for 50 ml rotor bores	1 Pack × 8 Pcs.	E1450-9015
	Storage box 5.0ml tubes, blue	3 Pcs. (1 Pack × 3 Pcs.)	E2325-5021

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk