



Standard screw cap for 5.0 ml tubes, white



Product Highlights

- 3/4 turn leakproof screw cap
- For use with 5.0 ml screw cap centrifuge tubes
- RNase, DNase, DNA & Pyrogen free
- Temperature use: -196°C to 100°C
- Centrifugation rate: 25,000 x g
- Cap center with "membrane" area for easy puncture with a needle
- Made of high density polyethylene
- Caps are not autoclavable

Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk




General Data

Art. No.	E1450-5060
Pack Size	200 Pcs. (2 Bags × 100 Pcs.)
Diameter	19.9 mm
Sterile	No
Autoclavable	No
Cap style	Standard screw cap
Material	Polyethylene
DNA free	Yes
DNase free	Yes
Pyrogen free	Yes
RNase free	Yes
Color	White
g-Force max.	25.000 x g

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk




All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 Standard screw cap for 5.0 ml tubes, white Cap style: Standard screw cap  Create PDF Data Sheet	200 Pcs. (2 Bags × 100 Pcs.)	E1450-5060
 Tethered screw cap for 5.0 ml tubes, white Cap style: Tethered screw cap  Create PDF Data Sheet	200 Pcs. (2 Bags × 100 Pcs.)	E1450-5080

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



Accessories

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	Rack 5.0ml tubes, blue	5 Pcs. (1 Pack × 5 Pcs.)	E2325-5011
	Universal centrifuge adapter, for 5.0 ml tubes, for 15 ml rotor bores	1 Pack × 8 Pcs.	E1450-9011
	Universal centrifuge adapter, for 5.0 ml tubes, for 50 ml rotor bores	1 Pack × 8 Pcs.	E1450-9015
	Storage box 5.0ml tubes, blue	3 Pcs. (1 Pack × 3 Pcs.)	E2325-5021

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk