



## Chest Freezer Racks for Standard Storage Boxes



### Product Highlights

- Racks for use in chest freezers have a locking rod to keep boxes secure and a handle on top
- For use with all standard storage boxes up to 136 x 136 x 53 mm
- Laser cutting technology results in high quality units with no sharp edges
- Made from 85% stainless steel
- Made in a factory powered 100% by sustainable energy
- Article No. etched onto every rack for easy reordering

## Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data

Art. No.	See variations
Rack Type	Chest Freezer Rack
Material	Stainless Steel

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## More informations about Chest Freezer Racks for Standard Storage Boxes











Chest [freezer](#) racks for use with all standard [storage boxes](#) up to 136 x 136 x 53 mm. Keeping your [freezer](#) tidy and optimising storage space is an efficient use of energy. Racks for chest freezers feature a locking rod to keep your [storage boxes](#) secure and a handle on top for easy removal from the [freezer](#)

The following [storage boxes](#) are suitable for use with these racks:

- StarStore 81, StarStore 49, StarStore
- The 50 mm high cardboard storage boxes
- 81-place cryogenic storage boxes, 100-place cryogenic storage boxes
- StarRack 64



## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">9-High Rack for Chest Freezer</a> <b>Places:</b> 9 (9 high)  <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Rack × 1 Piece)	E6010-0009
 <a href="#">10-High Rack for Chest Freezer</a> <b>Places:</b> 10 (10 high)  <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Rack × 1 Piece)	E6010-0010
 <a href="#">11-High Rack for Chest Freezer</a> <b>Places:</b> 11 (11 high)  <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Rack × 1 Piece)	E6010-0011
 <a href="#">12-High Rack for Chest Freezer</a> <b>Places:</b> 12 (12 high)  <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Rack × 1 Piece)	E6010-0012
 <a href="#">13-High Rack for Chest Freezer</a> <b>Places:</b> 13 (13 high)  <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Rack × 1 Piece)	E6010-0013

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)