



## Screw Cap Microtubes with Natural Tethered Caps (Sterile)



### Product Highlights

- Routinely tested to 20,000 x g
- Sterile and certified RNase, DNase, DNA, PCR Inhibitor and Endotoxin free.
- Tubes are supplied with pre-assembled natural tethered caps
- Tethered cap design ensures that the cap is always attached
- Tube made of polypropylene (PP)

## Product Applications

Suitable for many applications from sample storage in the gas phase of liquid nitrogen to freezer, through boiling to autoclave

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data

|                    |                               |
|--------------------|-------------------------------|
| Art. No.           | See variations                |
| Sterile            | Yes                           |
| Autoclavable       | No                            |
| Pack Size          | 250 Tubes (5 Bags × 50 Tubes) |
| Cap style          | Natural tethered cap          |
| Material           | Polypropylene                 |
| DNA free           | Yes                           |
| DNase free         | Yes                           |
| Free of endotoxins | Yes                           |
| PCR inhibitor free | Yes                           |
| RNase free         | Yes                           |
| g-Force max.       | 20,000 x g                    |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)













## More informations about Screw Cap Microtubes with Natural Tethered Caps (Sterile)

### Assembled sterile tubes & caps

A range of sterile and pre-assembled microtubes, available with a range of different cap styles. Supplied in convenient, resealable mini bags which stand up on your bench. Suitable for many applications from sample storage in the gas phase of liquid nitrogen through [freezer](#) to boiling.



## All Variations

| PRODUCT NAME                                                                                                                                                                                                                                                                                                                                                                                      | PACKAGING SIZE                | ART. NO.   |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------|------------|
| <br><a href="#">0.5 ml Plain Skirted Tube, Natural Tethered Screw Cap (Sterile)</a><br><b>Volume:</b> 0.5 ml<br><b>Tube Base:</b> Skirted<br><b>Tube Style:</b> Plain<br><br> <a href="#">Create PDF Data Sheet</a>             | 250 Tubes (5 Bags × 50 Tubes) | E1405-2142 |
| <br><a href="#">1.5 ml Graduated Conical Tube, Natural Tethered Screw Cap (Sterile)</a><br><b>Volume:</b> 1.5 ml<br><b>Tube Base:</b> Conical<br><b>Tube Style:</b> Graduated<br><br> <a href="#">Create PDF Data Sheet</a>     | 250 Tubes (5 Bags × 50 Tubes) | E1415-2232 |
| <br><a href="#">1.5 ml Plain Skirted Tube, Natural Tethered Screw Cap (Sterile)</a><br><b>Volume:</b> 1.5 ml<br><b>Tube Base:</b> Skirted<br><b>Tube Style:</b> Plain<br><br> <a href="#">Create PDF Data Sheet</a>          | 250 Tubes (5 Bags × 50 Tubes) | E1415-2242 |
| <br><a href="#">2.0 ml Graduated Conical Tube, Natural Tethered Screw Cap (Sterile)</a><br><b>Volume:</b> 2.0 ml<br><b>Tube Base:</b> Conical<br><b>Tube Style:</b> Graduated<br><br> <a href="#">Create PDF Data Sheet</a> | 250 Tubes (5 Bags × 50 Tubes) | E1420-2332 |
| <br><a href="#">2.0 ml Graduated Skirted Tube, Natural Tethered Screw Cap (Sterile)</a><br><b>Volume:</b> 2.0 ml<br><b>Tube Base:</b> Skirted<br><b>Tube Style:</b> Graduated<br><br> <a href="#">Create PDF Data Sheet</a> | 250 Tubes (5 Bags × 50 Tubes) | E1420-2342 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)