



## Screw Cap Microtubes with White/Colour EasyGrip Caps (Sterile)



### Product Highlights

- Routinely tested to 20,000 x g
- Sterile and certified RNase, DNase, DNA, PCR Inhibitor and Endotoxin free.
- Tubes supplied ready assembled with natural EasyGrip Caps with white integral grip/seal.
- One piece construction EasyGrip cap ensures a leak-free seal
- Tube made of polypropylene (PP)

## Product Applications

Suitable for many applications from sample storage in the gas phase of liquid nitrogen to freezer, through boiling to autoclave

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data

Art. No.	See variations
Sterile	Yes
Autoclavable	No
Pack Size	250 Tubes (5 Bags × 50 Tubes)
Cap style	White/colour – EasyGrip caps
Material	Polyethylene
DNA free	Yes
DNase free	Yes
Free of endotoxins	Yes
PCR inhibitor free	Yes
RNase free	Yes
Color	Various
g-Force max.	20,000 x g

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## More informations about Screw Cap Microtubes with White/Colour EasyGrip Caps (Sterile)

### Assembled sterile tubes & caps











A range of sterile and pre-assembled microtubes, available with a range of different cap styles. Supplied in convenient, resealable mini bags which stand up on your bench. Suitable for many applications from [sample storage](#) in the gas phase of liquid nitrogen through [freezer](#) to boiling.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)















## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">0.5 ml Plain Skirted Tube, EasyGrip Screw Cap White/Blue (Sterile)</a> <b>Volume:</b> 0.5 ml <b>Tube Base:</b> Skirted <b>Color:</b> Blue   <a href="#">Create PDF Data Sheet</a>	250 Tubes (5 Bags × 50 Tubes)	E1405-0301
 <a href="#">0.5 ml Plain Skirted Tube with EasyGrip Screw Cap White/Red (Sterile)</a> <b>Volume:</b> 0.5 ml <b>Tube Base:</b> Skirted <b>Color:</b> Red   <a href="#">Create PDF Data Sheet</a>	250 Tubes (5 Bags × 50 Tubes)	E1405-0304
 <a href="#">0.5 ml Plain Skirted Tube with EasyGrip Screw Cap White/Yellow (Sterile)</a> <b>Volume:</b> 0.5 ml <b>Tube Base:</b> Skirted <b>Color:</b> Yellow   <a href="#">Create PDF Data Sheet</a>	250 Tubes (5 Bags × 50 Tubes)	E1405-0306
 <a href="#">0.5 ml Plain Skirted Tube with EasyGrip Screw Cap White/Green (Sterile)</a> <b>Volume:</b> 0.5 ml <b>Tube Base:</b> Skirted <b>Color:</b> Green   <a href="#">Create PDF Data Sheet</a>	250 Tubes (5 Bags × 50 Tubes)	E1405-0302
 <a href="#">1.5 ml Graduated Conical Tube with EasyGrip Screw Cap White/Blue (Sterile)</a> <b>Volume:</b> 1.5 ml <b>Tube Base:</b> Conical <b>Color:</b> Blue   <a href="#">Create PDF Data Sheet</a>	250 Tubes (5 Bags × 50 Tubes)	E1415-0301

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)















PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">1.5ml Graduated Conical Tube with EasyGrip Screw Cap White/Green (Sterile)</a> <b>Volume:</b> 1.5 ml <b>Tube Base:</b> Conical <b>Color:</b> Green   <a href="#">Create PDF Data Sheet</a>	250 Tubes (5 Bags × 50 Tubes)	E1415-0302
 <a href="#">1.5ml Graduated Conical Tube with EasyGrip Screw Cap White/Red (Sterile)</a> <b>Volume:</b> 1.5 ml <b>Tube Base:</b> Conical <b>Color:</b> Red   <a href="#">Create PDF Data Sheet</a>	250 Tubes (5 Bags × 50 Tubes)	E1415-0304
 <a href="#">1.5ml Graduated Conical Tube with EasyGrip Screw Cap White/Yellow (Sterile)</a> <b>Volume:</b> 1.5 ml <b>Tube Base:</b> Conical <b>Color:</b> Yellow   <a href="#">Create PDF Data Sheet</a>	250 Tubes (5 Bags × 50 Tubes)	E1415-0306
 <a href="#">1.5ml Plain Skirted Tube with EasyGrip Screw Cap White/Blue (Sterile)</a> <b>Volume:</b> 1.5 ml <b>Tube Base:</b> Skirted <b>Color:</b> Blue   <a href="#">Create PDF Data Sheet</a>	250 Tubes (5 Bags × 50 Tubes)	E1415-0311
 <a href="#">1.5ml Plain Skirted Tube with EasyGrip Screw Cap White/Green (Sterile)</a> <b>Volume:</b> 1.5 ml <b>Tube Base:</b> Skirted <b>Color:</b> Green   <a href="#">Create PDF Data Sheet</a>	250 Tubes (5 Bags × 50 Tubes)	E1415-0312
 <a href="#">1.5ml Plain Skirted Tube with EasyGrip Screw Cap White/Red (Sterile)</a> <b>Volume:</b> 1.5 ml <b>Tube Base:</b> Skirted <b>Color:</b> Red   <a href="#">Create PDF Data Sheet</a>	250 Tubes (5 Bags × 50 Tubes)	E1415-0314

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)









PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <p><a href="#">1.5ml Plain Skirted Tube with EasyGrip Screw Cap White/Yellow (Sterile)</a> <b>Volume:</b> 1.5 ml <b>Tube Base:</b> Skirted <b>Color:</b> Yellow</p> <p> <a href="#">Create PDF Data Sheet</a></p>	250 Tubes (5 Bags × 50 Tubes)	E1415-0316
 <p><a href="#">2.0ml Graduated Conical Tube with EasyGrip Screw Cap White/Blue (Sterile)</a> <b>Volume:</b> 2.0 ml <b>Tube Base:</b> Conical <b>Color:</b> Blue</p> <p> <a href="#">Create PDF Data Sheet</a></p>	250 Tubes (5 Bags × 50 Tubes)	E1420-0301
 <p><a href="#">2.0ml Graduated Conical Tube with EasyGrip Screw Cap White/Green (Sterile)</a> <b>Volume:</b> 2.0 ml <b>Tube Base:</b> Conical <b>Color:</b> Green</p> <p> <a href="#">Create PDF Data Sheet</a></p>	250 Tubes (5 Bags × 50 Tubes)	E1420-0302
 <p><a href="#">2.0ml Graduated Conical Tube with EasyGrip Screw Cap White/Red (Sterile)</a> <b>Volume:</b> 2.0 ml <b>Tube Base:</b> Conical <b>Color:</b> Red</p> <p> <a href="#">Create PDF Data Sheet</a></p>	250 Tubes (5 Bags × 50 Tubes)	E1420-0304
 <p><a href="#">2.0ml Graduated Conical Tube with EasyGrip Screw Cap White/Yellow (Sterile)</a> <b>Volume:</b> 2.0 ml <b>Tube Base:</b> Conical <b>Color:</b> Yellow</p> <p> <a href="#">Create PDF Data Sheet</a></p>	250 Tubes (5 Bags × 50 Tubes)	E1420-0306
 <p><a href="#">2.0ml Graduated Skirted Tube with EasyGrip Screw Cap White/Blue (Sterile)</a> <b>Volume:</b> 2.0 ml <b>Tube Base:</b> Skirted <b>Color:</b> Blue</p> <p> <a href="#">Create PDF Data Sheet</a></p>	250 Tubes (5 Bags × 50 Tubes)	E1420-0311

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">2.0ml Graduated Skirted Tube with EasyGrip Screw Cap White/Green (Sterile)</a> <b>Volume:</b> 2.0 ml <b>Tube Base:</b> Skirted <b>Color:</b> Green   <a href="#">Create PDF Data Sheet</a>	250 Tubes (5 Bags × 50 Tubes)	E1420-0312
 <a href="#">2.0ml Graduated Skirted Tube with EasyGrip Screw Cap White/Red (Sterile)</a> <b>Volume:</b> 2.0 ml <b>Tube Base:</b> Skirted <b>Color:</b> Red   <a href="#">Create PDF Data Sheet</a>	250 Tubes (5 Bags × 50 Tubes)	E1420-0314
 <a href="#">2.0ml Graduated Skirted Tube with EasyGrip Screw Cap White/Yellow (Sterile)</a> <b>Volume:</b> 2.0 ml <b>Tube Base:</b> Skirted <b>Color:</b> Yellow   <a href="#">Create PDF Data Sheet</a>	250 Tubes (5 Bags × 50 Tubes)	E1420-0316

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)