



8-Strip PCR Tubes



Product Highlights

- Manufactured from medical grade polypropylene
- Tubes can easily be split into smaller sections by snapping the strips apart
- Ultra thin walls ensure efficient thermal transfer and maximum yield
- PCR inhibitor free
- Certified RNase, DNase, DNA and Pyrogen-free
- Supplied in resealable bags

Product Applications

PCR

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



General Data

| | |
|-----------------------------------|----------------|
| Art. No. | See variations |
| Sterile | No |
| Autoclavable | Yes |
| Material | Polypropylene |
| DNA free | Yes |
| DNase free | Yes |
| Human DNA free | Yes |
| Pyrogen free | Yes |
| RNase free | Yes |
| Color | Various |
| Suitable for PCR | yes |
| Suitable for Real-Time PCR (qPCR) | yes |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



More informations about 8-Strip PCR Tubes











8-strip PCR tubes. A range of different types of PCR strip caps are available separately, so you can mix and match tube and cap style to suit your needs.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk





All Variations

| PRODUCT NAME | PACKAGING SIZE | ART. NO. |
|--|---|--------------|
|  0.2 ml 8-Strip PCR Tubes, Natural Volume: 0.2 ml Color: Natural Pack Size: 125 Tube Strips (1 Bag × 125 Tube Strips)  Create PDF Data Sheet | 125 Tube Strips (1 Bag × 125 Tube Strips) | I1402-3500 |
|  0.2 ml 8-Strip PCR Tubes (Case Size), natural Volume: 0.2 ml Color: Natural Pack Size: 1,250 Tube Strips (10 Bags × 125 Tube Strips)  Create PDF Data Sheet | 1,250 Tube Strips (10 Bags × 125 Tube Strips) | I1402-3500-C |
|  0.2 ml 8-Strip PCR Tubes, Mixed Volume: 0.2 ml Color: Mixed Pack Size: 125 Tube Strips (5 Bags × 25 Tube Strips)  Create PDF Data Sheet | 125 Tube Strips (5 Bags × 25 Tube Strips) | I1402-3508 |
|  0.2 ml 8-Strip PCR Tubes, White Volume: 0.2 ml Color: White Pack Size: 125 Tube Strips (1 Bag × 125 Tube Strips)  Create PDF Data Sheet | 125 Tube Strips (1 Bag × 125 Tube Strips) | I1402-3509 |
|  0.1 ml 8-Strip PCR Tubes, Low-Profile, Natural Volume: 0.1 ml Color: Natural Pack Size: 125 Tube Strips (1 Bag × 125 Tube Strips)  Create PDF Data Sheet | 125 Tube Strips (1 Bag × 125 Tube Strips) | I1402-3600 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



| PRODUCT NAME | PACKAGING SIZE | ART. NO. |
|--|---|------------|
|  <p>0.1 ml 8-Strip PCR Tubes, Low-Profile, white Volume: 0.1 ml Color: White Pack Size: 125 Tube Strips (1 Bag × 125 Tube Strips)</p> <p> Create PDF Data Sheet</p> | 125 Tube Strips (1 Bag × 125 Tube Strips) | I1402-3609 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk