



## Cryogenic Lock Racks



### Product Highlights

- 50-place polypropylene rack in a 5x10 array for skirted cryogenic vials
- Skirted vials lock into the rack for easy one-handed cap removal
- Anti-skid rubber feet
- Temperature use: -90 °C to 121 °C
- Autoclavable
- Dimensions (WxDxH): 200 x 100 x 25 mm

## Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data

Art. No.	See variations
Dimensions (W × D × H)	200 mm × 100 mm × 25 mm
Dimensions (W × D)	200 mm × 100 mm
Autoclavable	Yes
Pack Size	4 Pcs. (4 Pcs.)
Material	Polypropylene
Color	Various
Places	50

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## More informations about Cryogenic Lock Racks







50-place polypropylene rack in three different colors for skirted cryogenic vials. The vials lock into the rack for easy one-handed cap removal.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">50 Place Lock-Rack, Blue</a> Color: Blue  <a href="#">Create PDF Data Sheet</a>	4 Pcs. (1 Pack × 4 Pcs.)	E2350-4001
 <a href="#">50 Place Lock-Rack, Red</a> Color: Red  <a href="#">Create PDF Data Sheet</a>	4 Pcs. (1 Pack × 4 Pcs.)	E2350-4004
 <a href="#">50 Place Lock-Rack, Violet</a> Color: Violet  <a href="#">Create PDF Data Sheet</a>	4 Pcs. (1 Pack × 4 Pcs.)	E2350-4005

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)