




## 12-Strip PCR Tubes



### Product Highlights

- Manufactured from medical grade polypropylene
- Tubes can easily be split into smaller sections by snapping the strips apart
- Ultra thin walls ensure efficient thermal transfer and maximum yield
- PCR inhibitor free
- Certified RNase, DNase, DNA and Pyrogen-free
- Supplied in resealable bags

## Product Applications

PCR

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data

|                                   |                |
|-----------------------------------|----------------|
| Art. No.                          | See variations |
| Sterile                           | No             |
| Autoclavable                      | Yes            |
| Volume                            | 0.2 ml         |
| Material                          | Polypropylene  |
| DNA free                          | Yes            |
| DNase free                        | Yes            |
| Human DNA free                    | Yes            |
| PCR inhibitor free                | Yes            |
| Pyrogen free                      | Yes            |
| RNase free                        | Yes            |
| Color                             | Various        |
| Suitable for PCR                  | yes            |
| Suitable for Real-Time PCR (qPCR) | yes            |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## More informations about 12-Strip PCR Tubes

8- and 12-strip PCR tubes. A range of different types of PCR strip caps are available separately, so you can mix and match tube and cap style to suit your needs.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## All Variations

| PRODUCT NAME   | PACKAGING SIZE                             | ART. NO.     |
|--|--|--------------|
| <br><a href="#">0.2 ml 12-Strip PCR Tubes, Natural</a><br><b>Color:</b> Natural<br><b>Pack Size:</b> 80 Tube Strips<br>(1 Bag × 80 Tube Strips)<br> <a href="#">Create PDF Data Sheet</a>                | 80 Tube Strips (1 Bag × 80 Tube Strips)    | I1402-1200   |
| <br><a href="#">0.2 ml 12-Strip PCR Tubes, Natural (Case Size)</a><br><b>Color:</b> Natural<br><b>Pack Size:</b> 800 Tube Strips<br>(10 Bags × 80 Tube Strips)<br> <a href="#">Create PDF Data Sheet</a> | 800 Tube Strips (10 Bags × 80 Tube Strips) | I1402-1200-C |
| <br><a href="#">0.2 ml 12-Strip PCR Tubes, Mixed</a><br><b>Color:</b> Mixed<br><b>Pack Size:</b> 80 Tube Strips<br>(5 Bags × 16 Tube Strips)<br> <a href="#">Create PDF Data Sheet</a>                | 80 Tube Strips (5 Bags × 16 Tube Strips)   | I1402-1208   |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)