




8-Strip PCR Caps



Product Highlights

- Manufactured from medical grade polypropylene
- 8-Strip Caps can be easily split into smaller sections by snapping apart
- Certified RNase, DNase, DNA and Pyrogen-free
- Supplied in resealable bags

Product Applications

PCR

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



General Data

Art. No.	See variations
Sterile	No
Autoclavable	Yes
Material	Polypropylene
DNA free	Yes
DNase free	Yes
Free of endotoxins	Yes
Human DNA free	Yes
PCR inhibitor free	Yes
Pyrogen free	Yes
RNase free	Yes
Color	Various
Suitable for PCR	yes
For use with (SL)	PCR Tubes & Plates

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



More informations about 8-Strip PCR Caps





- › PCR tube and cap strips made from medical grade polypropylene
- › Strips can be easily split into smaller sections by snapping apart
- › Extra clear flat cap strips suitable for qPCR
- › Supplied in resealable bags
- › Certified RNase, DNase, DNA and Pyrogen-free
- › PCR Inhibitor free

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk





All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <p>8-Strip PCR Caps, Domed, Natural Cap style: Domed Color: Natural Pack Size: 125 Pcs. (1 Bag × 125 Pcs.)</p> <p> Create PDF Data Sheet</p>	125 Pcs. (1 Bag × 125 Pcs.)	I1400-0800
 <p>8-Strip PCR Caps, Domed, Natural (Case Size) Cap style: Domed Color: Natural Pack Size: 1,250 Pcs. (10 Bags × 125 Pcs.)</p> <p> Create PDF Data Sheet</p>	1,250 Pcs. (10 Bags × 125 Pcs.)	I1400-0800-C
 <p>8-Strip PCR Caps, Domed, Mixed Cap style: Domed Color: Mixed Pack Size: 125 Pcs. (5 Bags × 25 Pcs.)</p> <p> Create PDF Data Sheet</p>	125 Pcs. (5 Bags × 25 Pcs.)	I1400-0808
 <p>8-Strip PCR Caps, Flat (Xtra-Clear), Natural Cap style: Flat (Xtra-Clear) Color: Natural Pack Size: 125 Pcs. (1 Bag × 125 Pcs.)</p> <p> Create PDF Data Sheet</p>	125 Pcs. (1 Bag × 125 Pcs.)	I1400-0900
 <p>8-Strip PCR Caps, Flat (Xtra-Clear), Natural (Case Size) Cap style: Flat (Xtra-Clear) Color: Natural Pack Size: 1,250 Pcs. (10 Bags × 125 Pcs.)</p> <p> Create PDF Data Sheet</p>	1,250 Pcs. (10 Bags × 125 Pcs.)	I1400-0900-C
 <p>8-Strip Caps, Domed, Perfect Seal, Natural Cap style: Domed, perfect seal Color: Natural Pack Size: 125 Pcs. (1 Bag × 125 Pcs.)</p> <p> Create PDF Data Sheet</p>	125 Pcs. (1 Bag × 125 Pcs.)	I1400-0810

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



PRODUCT NAME	PACKAGING SIZE	ART. NO.
 8-Strip Caps, Domed, Perfect Seal, Mixed Cap style: Domed, perfect seal Color: Mixed Pack Size: 125 Pcs. (5 Bags × 25 Pcs.)  Create PDF Data Sheet	125 Pcs. (5 Bags × 25 Pcs.)	I1400-0818

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk