



## Semi-Skirted PCR Plates with Raised Rim



### Product Highlights

- ABI FAST®-Style plate
- Certified RNase, DNase, DNA and PCR Inhibitor free
- Tested Pyrogen free

## Product Applications

PCR

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## General Data

Art. No.	See variations
Sterile	No
Autoclavable	No
Number of wells	96 wells
Material	Polypropylene
DNA free	Yes
DNase free	Yes
Free of endotoxins	Yes
Free of heavy metals	Yes
PCR inhibitor free	Yes
RNase free	Yes
Frame Design	Semi-skirted
Suitable for PCR	yes
Suitable for Real-Time PCR (qPCR)	yes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)













## More informations about Semi-Skirted PCR Plates with Raised Rim

Starlab's semi-skirted [PCR plates](#) with raised rim are made from medical grade polypropylene to provide the quality assurance you need. More rigid than non-skirted plates, they also offer better handling convenience. The plates have clearly marked wells to provide easy identification and sample tracking, as well as compatibility with a wide range of ABI cyclers and sequencers.

Supplied in convenient, resealable bags to minimise contamination risk. Certified free of RNase, DNase, DNA and PCR Inhibitors. Tested Pyrogen free.





## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">96-Well PCR Plate, Semi-Skirted, Raised Rim, 0.1ml Low Profile (ABI FAST®-Style), natural</a> <b>Profile:</b> Low <b>Color:</b> Natural <b>Pack Size:</b> 10 Plates (1 Box × 10 Plates)  <a href="#">Create PDF Data Sheet</a>	10 Plates (1 Box × 10 Plates)	E1403-7700
 <a href="#">96-Well PCR Plate, Semi-Skirted, Raised Rim, 0.1ml Low Profile (ABI FAST®-Style), natural (Case Size)</a> <b>Profile:</b> Low <b>Color:</b> Natural <b>Pack Size:</b> 100 Plates (10 Boxes × 10 Plates)  <a href="#">Create PDF Data Sheet</a>	100 Plates (10 Boxes × 10 Plates)	E1403-7700-C
 <a href="#">96-Well PCR Plate, Semi-Skirted, Raised Rim, 0.1ml Low Profile (ABI FAST®-Style), white</a> <b>Profile:</b> Low <b>Color:</b> White <b>Pack Size:</b> 10 Plates (1 Box × 10 Plates)  <a href="#">Create PDF Data Sheet</a>	10 Plates (1 Box × 10 Plates)	E1403-7709
 <a href="#">96-Well PCR Plate, Semi-Skirted, Raised Rim, 0.1ml Low Profile (ABI FAST®-Style), white (Case Size)</a> <b>Profile:</b> Low <b>Color:</b> White <b>Pack Size:</b> 100 Plates (10 Boxes × 10 Plates)  <a href="#">Create PDF Data Sheet</a>	100 Plates (10 Boxes × 10 Plates)	E1403-7709-C
 <a href="#">96-Well PCR Plate, Semi-Skirted, Raised Rim, 0.2ml (ABI-Style), natural</a> <b>Profile:</b> Standard height <b>Color:</b> Natural <b>Pack Size:</b> 10 Plates (1 Box × 10 Plates)  <a href="#">Create PDF Data Sheet</a>	10 Plates (1 Box × 10 Plates)	E1403-8200

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <p><a href="#">96-Well PCR Plate, Semi-Skirted, Raised Rim, 0.2ml (ABI-Style), natural (Case Size)</a> <b>Profile:</b> Standard height <b>Color:</b> Natural <b>Pack Size:</b> 100 Plates (10 Boxes × 10 Plates)</p> <p> <a href="#">Create PDF Data Sheet</a></p>	100 Plates (10 Boxes × 10 Plates)	E1403-8200-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)