



## 50 ml Centrifuge Tubes



### Product Highlights

- Clear, sterile polypropylene tubes with green HDPE screw caps
- White writing surface and clear, printed graduations
- Perfect for centrifugation of cells: USP Class VI, non-cytotoxic and sterile (SAL  $10^{-6}$ ) material
- Certified RNase, DNase, DNA, heavy-metal and Pyrogen free

## Product Applications

Centrifugation  
Sample preparation  
Preparation of master mixes and buffers  
Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data

Art. No.	See variations
Dimensions (W × H)	34.9 mm × 115.5 mm
Sterile	Yes
Autoclavable	Yes
Volume	50 ml
Material	Polypropylene
DNA free	Yes
DNase free	Yes
Free of endotoxins	Yes
Free of heavy metals	Yes
Human DNA free	Yes
Pyrogen free	Yes
RNase free	Yes
Color	Clear
g-Force max.	12,000 x g
Non-cytotoxic	Yes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)









## More informations about 50 ml Centrifuge Tubes

STARLAB's new [centrifuge tubes](#) are designed for the centrifugation of cells, but can also be used for standard centrifugation, preparation and storage applications. USP Class VI polymers and certified non-toxicity ensure the viability and integrity of cells. The sterility of assurance level (SAL) of  $10^{-6}$  fulfills highest demands on sterility, meeting EN 556-1 standard (Sterilisation of medical devices). Racked tubes are delivered in recyclable cardboard racks.

**Maximum centrifugation rates:** Conical Tubes 15 & 50 ml: 12,000 x g and Skirted Tube 50 ml: 6,000 x g



## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">50 ml Centrifuge Tube, Conical, Loose (Sterile)</a> <b>Tube Base:</b> Conical <b>Packaging Unit:</b> Bags   <a href="#">Create PDF Data Sheet</a>	500 Tubes (20 Bags × 25 Tubes)	E1450-0100
 <a href="#">50 ml Centrifuge Tube, Conical, Racked (Sterile)</a> <b>Tube Base:</b> Conical <b>Packaging Unit:</b> Racks   <a href="#">Create PDF Data Sheet</a>	500 Tubes (20 Racks × 25 Tubes)	E1450-0700
 <a href="#">50 ml Centrifuge Tube, Skirted, Loose (Sterile)</a> <b>Tube Base:</b> Skirted <b>Packaging Unit:</b> Bags   <a href="#">Create PDF Data Sheet</a>	500 Tubes (20 Bags × 25 Tubes)	E1450-0400

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)