



ErgoOne® E 8-Channel Electronic Pipette, 50 – 1250µl



Product Highlights

- Easy to use, intuitive menu structure
- Lightweight, ergonomic design
- Adjustable finger rest for individual comfort
- Significantly reduced tip loading and ejection forces
- Universal tip fit
- Reliable
- Autoclavable pipette shaft
- 12-month warranty
- Includes two free racks of matching TipOne® tips
- Includes shelf holder G9999-1009
- supplied with battery, mains adapter, shelf holder (G9999-1009), reagent reservoir, 1 set of O-rings, 2 racks with tips and silicone grease

Product Applications

Pipetting

Reverse Pipetting of "tricky" liquids (e. g. with high viscosity)

Mixing

Dispense

Electrophoresis (excl. 1250 µl models)

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



General Data

Art. No.	G9008-1250
Pack Size	1 Piece (1 Box × 1 Piece)
Color code	Blue
Autoclavable	Autoclavable pipette shaft
Volume range	50 - 1,250 µl
Volume	1250 µl
Volume Setting	Variable
Channel type	8-channel
Operating mode	Electronic
Pipetting type	Air-cushion
Operation	Battery operated piston stroke / manual tip ejection

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



More informations about ErgoOne® E 8-Channel Electronic Pipette, 50 – 1250µl

Reliable and ergonomic!













ErgoOne® E multi-channel pipettes provide reliable pipetting every day. Ergonomic and lightweight, ErgoOne® E pipettes are perfect for comfortable pipetting during long periods of use.

Five operating modes provide standard and reverse pipetting, mixing and dispensing, plus a unique “Electrophoresis” mode to help you load gels easily. The long lasting battery can also be charged during use so you are never without your pipette.

Please note that GEL mode is not available for ErgoOne® E 1000 µl, 1250 µl and 5000 µl











All Variations

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	ErgoOne® E 8-Channel Electronic Pipette, 0.5 – 10µl Channel type: 8-channel Volume range: 0.5 - 10 µl  Create PDF Data Sheet	1 Piece (1 Box × 1 Piece)	G9008-0010
	ErgoOne® E 8-Channel Electronic Pipette, 1.0 – 20µl Channel type: 8-channel Volume range: 1.0 - 20 µl  Create PDF Data Sheet	1 Piece (1 Box × 1 Piece)	G9008-0020
	ErgoOne® E 8-Channel Electronic Pipette, 5.0 – 100µl Channel type: 8-channel Volume range: 5.0 - 100 µl  Create PDF Data Sheet	1 Piece (1 Box × 1 Piece)	G9008-0100
	ErgoOne® E 8-Channel Electronic Pipette, 15 – 300µl Channel type: 8-channel Volume range: 15 - 300 µl  Create PDF Data Sheet	1 Piece (1 Box × 1 Piece)	G9008-0300
	ErgoOne® E 8-Channel Electronic Pipette, 50 – 1250µl Channel type: 8-channel Volume range: 50 - 1,250 µl  Create PDF Data Sheet	1 Piece (1 Box × 1 Piece)	G9008-1250
	ErgoOne® E 12-Channel Pipette Electronic, 5.0 – 100µl Channel type: 12-channel Volume range: 5.0 - 100 µl  Create PDF Data Sheet	1 Piece (1 Box × 1 Piece)	G9012-0100

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.


Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk














PRODUCT NAME	PACKAGING SIZE	ART. NO.
 ErgoOne® E 12-Channel Electronic Pipette, 0.5 – 10µl Channel type: 12-channel Volume range: 0.5 - 10 µl  Create PDF Data Sheet	1 Piece (1 Box × 1 Piece)	G9012-0010
 ErgoOne® E 12-Channel Electronic Pipette, 1.0 – 20µl Channel type: 12-channel Volume range: 1.0 - 20 µl  Create PDF Data Sheet	1 Piece (1 Box × 1 Piece)	G9012-0020
 ErgoOne® E 12-Channel Electronic Pipette, 15 – 300µl Channel type: 12-channel Volume range: 15 - 300 µl  Create PDF Data Sheet	1 Piece (1 Box × 1 Piece)	G9012-0300
 ErgoOne® E 12-Channel Electronic Pipette, 50 – 1250µl Channel type: 12-channel Volume range: 50 - 1,250 µl  Create PDF Data Sheet	1 Piece (1 Box × 1 Piece)	G9012-1250

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



Accessories

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 Shelf Holder for one ErgoOne® E pipette	1 Piece (1 Box × 1 Piece)	G9999-1009
 Wall Mount for one ErgoOne® E pipette	1 Piece (1 Box × 1 Piece)	G9999-1010
 Universal stand for one ErgoOne® E pipette - all volumes	1 Piece (1 Box × 1 Piece)	G9999-1011
 O-ring ErgoOne® E 8/12-c 1250 µl	1 Piece (1 Pack × 1 Piece)	P7121-1201
 Piston-shaft unit ErgoOne® E 1250 µl	1 Piece (1 Pack × 1 Piece)	P7141-1211
 Piston grease for ErgoOne® E 8- and 12-channel pipettes 1250 µl	1 Piece (1 Box × 1 Piece)	G9999-1012
 Fluorstatic Grease for Guiding Axis (for Multi-Channel Pipettes)	1 Piece (1 Pack × 1 Piece)	G9999-1008
 Silicone Oil for Piston (for 5000 µl Single-Channel and all Multi-Channel Pipettes)	1 Piece (1 Pack × 1 Piece)	G9999-1007
 Replacement Battery for ErgoOne® E	1 Piece (1 Box × 1 Piece)	G9999-1005
 AC Adapter for Europe/UK (230V/50 Hz)	1 Piece (1 Box × 1 Piece)	G9999-1000
 Individual Stand (for all ErgoOne® E Multi-Channel Pipettes up to 300 µl)	1 Piece (1 Box × 1 Piece)	G9999-1004

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk