



## 5 ml Gilson-Style Graduated Macro Tip (sterile)



### Product Highlights

- From 1200 µl tips to 10 ml tips
- Made from polypropylene
- For use with: ErgoOne®, Gilson®, Nichiryo®, StarPet, Volac®
- Different pack types available
- Certified
- DNase free, RNase free, DNA free, Pyrogene free
- DNase Testing is by gel electrophoresis. Test sensitivity 10<sup>-7</sup> Kunitz units/µl.
- RNase Testing is by gel electrophoresis. Test sensitivity 10<sup>-9</sup> Kunitz units/µl.
- DNA Testing is by gel electrophoresis. Test sensitivity 30 pg.
- Pyrogens Non-sterile product is tested pyrogen-free. Sterile product is tested and certified
- free of pyrogens/endotoxins by the LAL (Limulus Amebocyte Lysate) Gel Clot method. Test sensitivity 0.06 EU/ml (Endotoxin Units)
- PCR consumables are PCR inhibitor free

## Product Applications

Pipetting of large volumes of liquids

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## General Data

Art. No.	I1050-0710
Pack Size	500 Tips (10 Racks × 50 Tips)
Sterile	Yes
Autoclavable	Yes
Filter tip	Non-filter
Volume	5 ml
Material	Polypropylene
DNA free	Yes
DNase free	Yes
Pyrogen free	Yes
RNase free	Yes
Graduations	Yes
Design Tips	Gilson-Style
For use with	<b>Abimed</b> : Discovery Comfort, LABMATE Optima <b>Eppendorf</b> : Reference 2 <b>Gilson</b> : Pipetman Concept, Pipetman L, Pipetman P, Pipetman Ultra <b>HTL</b> : Discovery Comfort, Labmate Pro, Optipette <b>Nichiryo</b> : 31N, NPLO, NPX <b>Rainin</b> : Classic, EDP1, EDP-3 Plus, Pipet-Lite, Pipet-Plus <b>Socorex</b> : Acura electro 935 <b>STARLAB</b> : ErgoOne E, StarPet

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## More informations about 5 ml Gilson-Style Graduated Macro Tip (sterile)













Large Volume Tips

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)










## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">5 ml Biohit-Style Tip (non-sterile)</a> <b>Design Tips:</b> Biohit-Style <b>Sterile:</b> No <b>Pack Size:</b> 250 Tips (1 Bag × 250 Tips)  <a href="#">Create PDF Data Sheet</a>	250 Tips (1 Bag × 250 Tips)	E1009-7000
 <a href="#">5 ml Biohit-Style-XL Tip (non-sterile)</a> <b>Design Tips:</b> Biohit-Style-XL <b>Sterile:</b> No <b>Pack Size:</b> 250 Tips (1 Bag × 250 Tips)  <a href="#">Create PDF Data Sheet</a>	250 Tips (1 Bag × 250 Tips)	E1009-7001
 <a href="#">5 ml Finn-Style Tip (non-sterile)</a> <b>Design Tips:</b> Finn-Style <b>Sterile:</b> No <b>Pack Size:</b> 250 Tips (1 Bag × 250 Tips)  <a href="#">Create PDF Data Sheet</a>	250 Tips (1 Bag × 250 Tips)	I1009-5000
 <a href="#">5 ml Gilson-Style Graduated Macro Tip (non-sterile)</a> <b>Design Tips:</b> Gilson-Style <b>Sterile:</b> No <b>Pack Size:</b> 250 Tips (1 Bag × 250 Tips)  <a href="#">Create PDF Data Sheet</a>	250 Tips (1 Bag × 250 Tips)	I1050-0000
 <a href="#">5 ml Gilson-Style Graduated Macro Tip (non-sterile)</a> <b>Design Tips:</b> Gilson-Style <b>Sterile:</b> No <b>Pack Size:</b> 500 Tips (10 Racks × 50 Tips)  <a href="#">Create PDF Data Sheet</a>	500 Tips (10 Racks × 50 Tips)	I1050-0700
 <a href="#">5 ml Gilson-Style Graduated Macro Tip (sterile)</a> <b>Design Tips:</b> Gilson-Style <b>Sterile:</b> Yes <b>Pack Size:</b> 500 Tips (10 Racks × 50 Tips)  <a href="#">Create PDF Data Sheet</a>	500 Tips (10 Racks × 50 Tips)	I1050-0710

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">5 ml Eppendorf-Style Graduated Macro Tip (non-sterile)</a> <b>Design Tips:</b> Eppendorf-Style <b>Sterile:</b> No <b>Pack Size:</b> 250 Tips (1 Bag × 250 Tips)  <a href="#">Create PDF Data Sheet</a>	250 Tips (1 Bag × 250 Tips)	I1053-0000
 <a href="#">5 ml Eppendorf-Style Graduated Macro Tip (non-sterile)</a> <b>Design Tips:</b> Eppendorf-Style <b>Sterile:</b> No <b>Pack Size:</b> 500 Tips (10 Racks × 50 Tips)  <a href="#">Create PDF Data Sheet</a>	500 Tips (10 Racks × 50 Tips)	I1053-0700
 <a href="#">5 ml Eppendorf-Style Graduated Macro Tip (sterile)</a> <b>Design Tips:</b> Eppendorf-Style <b>Sterile:</b> Yes <b>Pack Size:</b> 500 Tips (10 Racks × 50 Tips)  <a href="#">Create PDF Data Sheet</a>	500 Tips (10 Racks × 50 Tips)	I1053-0710
 <a href="#">5.0 ml ErgoOne E/Finn-Style Pipette Tip (non-sterile), Loose</a> <b>Design Tips:</b> Finn-Style <b>Sterile:</b> No <b>Pack Size:</b> 200 Tips (1 Bag × 200 Tips)  <a href="#">Create PDF Data Sheet</a>	200 Tips (1 Bag × 200 Tips)	E1070-2596
 <a href="#">5.0 ml ErgoOne E/Finn-Style Pipette Tip (non-sterile), Loose</a> <b>Design Tips:</b> Finn-Style <b>Sterile:</b> No <b>Pack Size:</b> 1,000 Tips (5 Bags × 200 Tips)  <a href="#">Create PDF Data Sheet</a>	1,000 Tips (5 Bags × 200 Tips)	E1070-2601
 <a href="#">5 ml Gilson-Style Graduated Macro Tip (non-sterile)</a> <b>Design Tips:</b> Gilson-Style <b>Sterile:</b> No <b>Pack Size:</b> 2,500 Tips (10 Bag × 250 Tips)  <a href="#">Create PDF Data Sheet</a>	2,500 Tips (10 Bag × 250 Tips)	I1050-0000-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)