




## B-499 Nylon (25.40 x 12.70 x 11.18 mm), white



### Product Highlights

- Uses: Liquid nitrogen (-196°C) (not on glass), Freezer (-80°C), Autoclave (121°C)

## Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data

Art. No.	N0118-4999
Dimensions (W × H)	25.40 mm × 9.53 mm
External width	25.40 mm
Label Dimensions (mm) A	25.40 mm
Label Dimensions (mm) B	9.53 mm
Label Dimensions (mm) C	9.53 mm
Autoclave (121 °C)	Yes
Pack Size	350 Labels (1 Cartridge × 350 Labels)
Label Type	Free Standing
Label Style	2
Material	B-499 Nylon
Color	White
Compatible with 0.2 ml Tubes	no
Compatible with 0.5 / 0.65 ml Tubes	no
Compatible with 1.5/2.0 ml Tubes	yes
Compatible with 15 ml Tubes	no
Compatible with 1–8 ml Tubes	no
Compatible with 50 ml Tubes	no
Compatible with Plates	no
Compatible with Slides	no
Freezer (-80 °C)	Yes
Liquid Nitrogen (-196 °C)	Yes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



Slide Labelling

No

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## More informations about B-499 Nylon (25.40 x 12.70 x 11.18 mm), white

Label materials for laboratory applications such as cryogenic vials and slides. Ideal for general purpose labelling, the permanent adhesive and thicker construction makes this material durable yet easy to handle.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)