



## Chemgene MEDLAB Ready to Use Spray, Clear Unfragranced, 6 x 750 ml



### Product Highlights

- Cleans and disinfects to a consistently high standard
- Use against bacteria, mycobacteria, yeast, fungi, enveloped viruses, adenovirus and bacteriophage
- Pre-diluted solution (10%)
- Denatures RNA/DNA
- Suitable for use in laboratory, healthcare and forensic settings
- Not to be used to disinfect medical devices

## Product Applications

Disinfection  
Decontamination  
Cleaning  
Occupational health and safety

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data

Art. No.	SKU093A
Volume	750 ml
Pack Size	6 Trigger bottle (1 Box × 6 Trigger bottle)
Color	Clear
Shelf life	2 years
Fragrance (SL)	Unfragranced
Sustainability attributes	Spray bottle made made from 35 % post consumer recycled HDPE and is recyclable after use. Can also be reused. Recyclable shipper box made from 100 % recycled cardboard.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## More informations about Chemgene MEDLAB Ready to Use Spray, Clear Unfragranced, 6 x 750 ml

Chemgene MEDLAB multi-surface disinfectant delivers a new, unique formulation designed with laboratory, forensic and healthcare settings in mind.

The combination of active substances in Chemgene MEDLAB ensure that equipment and instruments can be routinely cleaned and disinfected to a consistently high standard, with less risk of antimicrobial resistance. Chemgene MEDLAB pre-diluted ready to use formula is effective from 60 seconds, and surfaces do not require rinsing following disinfection making the process quick and easy.

### PRODUCT FEATURES

- › Reliable pre-diluted solution (10%)
- › Denatures RNA/DNA
- › Effective against MDR pathogens
- › Non-corrosive
- › Biodegradable ingredients
- › Alcohol & chlorine free
- › No rinsing required

Chemgene MEDLAB Multi-Surface Disinfectants are suitable for the cleaning and disinfection of a variety of surfaces, instruments and equipment against bacteria, yeast, fungi, enveloped viruses, adenovirus and bacteriophage.

Use in the laboratory for protection from pathogenic microorganisms, in healthcare settings for protection against multi-drug resistant pathogens and forensic settings for protection from cross-contamination.

Chemgene MEDLAB multi-surface disinfectant is compliant with regulations in the UK and EU including the Biocidal Products Regulation (BPR), Classification of Labelling and Packaging (CLP) and Registration, Evaluation, Authorisation and Restriction of Chemicals (REACH) regulations.