



# Chemgene MEDLAB Reusable, Empty, Trigger Spray Bottle, 1 L



## Product Highlights

- Use for the compliant storage of diluted disinfectant
- Reusable, empty spray bottle with ergonomic trigger spray
- 1 L volume makes dilution of concentrate quick and easy
- Made from 100 % post-consumer HDPE

## Product Applications

- Disinfection
- Decontamination
- Cleaning
- Occupational health and safety

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data

Art. No.	SKU084A
Volume	1 l
Pack Size	5 Trigger bottle (1 Box × 5 Trigger bottle)
Solution	Empty, reusable bottle
Sustainability attributes	Reusable spray bottle made made from 100 % post consumer recycled HDPE, recyclable after use. Recyclable shipper box made from 100% recycled cardboard.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## More informations about Chemgene MEDLAB Reusable, Empty, Trigger Spray Bottle, 1 L

This reusable, empty Trigger Spray bottle holds 1 L of solution, making dilution of Chemgene MEDLAB from concentrate quick and easy.

### PRODUCT FEATURES

- > Reuse and refill with diluted concentrate
- > 1L volume makes dilution quick and easy
- > Ergonomic trigger spray
- > Made from 100 % post-consumer HDPE
- > Can be recycled at the end of its life

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)