



RepeatOne® S Dispenser



Product Highlights

- Up to 49 dispenses from the same sample
- Easy tip insertion – tip is simply inserted from below
- Increased chemical resistance through innovative plastic materials
- Conforms to EN ISO 8655 Piston-Operated Volumetric Apparatus
- Conforms to ASTM E1154-89 (2003) Standard Specification for Piston or Plunger Operated Volumetric Apparatus
- Slim, ergonomic and lightweight
- Convenient table of volume setting on the back
- Shelf mount included
- PFAs free
- Made in Germany, 1-year warranty

Product Applications

High-precision repetitive dispensing regardless of viscosity, vapour pressure and density of the liquid

Safe dispensing by avoidance of aerosol contaminations

Ideal for the dispensing of radioactive, infectious or toxic substances

Quick dispensing of long series

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



General Data

Art. No.	G4026-1500
Pack Size	1 Piece (1 Box × 1 Piece)
Sterile	No
Autoclavable	No
Volume range	2.0 µl - 5.0 ml
Volume	5.0 ml
Channel type	1-channel
Operating mode	Manual
Pipetting type	Positive displacement

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



More informations about RepeatOne® S Dispenser

Rapid dispensing!

RepeatOne® S Dispenser is a hand-held device for repetitive and rapid dispensing of volumes 2 µl to 5 ml.

Compatibility

RepeatOne® can be used with the following dispenser tips: Starlab TipOne® Repeat Dispenser Tips, Eppendorf, Brand, Gilson, HTL, Ritter.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



Accessories

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 Shelf Mount for RepeatOne® S Dispenser	1 Piece (1 Pack × 1 Piece)	G4026-1009

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk