



## Tubes for Microdilution Tube System



### Product Highlights

- For use with the reusable racks of the Microdilution Tube system
- Temperature use -80 °C to 121 °C
- Autoclavable
- Certified RNase, DNase, DNA and Pyrogen free

## Product Applications

Multichannel and automated pipetting

Preparation of serial dilutions

Mixing

Storage and harvesting of cells

Cell growth assays

DNA screening

Long-term storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data







Art. No.	See variations
Dimensions (W × H)	8.5 mm × 44.5 mm
Sterile	No
Autoclavable	Yes
Volume	1.2 ml
Material	Polypropylene
DNA free	Yes
DNase free	Yes
Pyrogen free	Yes
RNase free	Yes
Operating temperature	-80°C to 121°C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">1.2 ml Individual Reaction Tubes</a> <b>Tubes per strip:</b> individual  <a href="#">Create PDF Data Sheet</a>	1,000 Tubes (1 Bag × 1000 Tubes)	E1710-0000
 <a href="#">1.2 ml 8-Strip Reaction Tubes</a> <b>Tubes per strip:</b> 8  <a href="#">Create PDF Data Sheet</a>	125 Tubes (1 Bag × 125 Tubes)	E1720-8000
 <a href="#">1.2 ml 12-Strip Reaction Tubes</a> <b>Tubes per strip:</b> 12  <a href="#">Create PDF Data Sheet</a>	80 Tubes (1 Bag × 80 Tubes)	E1730-9000

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)