




Cardboard Storage Boxes for Centrifuge Tubes



Product Highlights

- Moisture-resistant coating
- 16-place box for 50 ml tubes
- 36-place box for 15 ml tubes
- Removable dividers

Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



General Data

| | |
|---------------------------|-------------------------|
| Art. No. | See variations |
| Dimensions (W × D × H) | 146 mm × 146 mm × 76 mm |
| Dimensions (W × D) | 146 mm × 146 mm |
| Max. tube length | 140 mm |
| Autoclavable | No |
| Material | Coated Cardboard |
| Color | White |
| Sustainability attributes | Recyclable |
| For use with (SL) | Centrifuge Tubes |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



More informations about Cardboard Storage Boxes for Centrifuge Tubes





White cardboard storage boxes with lift-off lid for the storage of 15 ml or 50 ml centrifuge tubes. Moisture-resistant coating. Removable dividers. Box dimensions (WxDxH): 146 x 146 x 76 mm.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



All Variations

| PRODUCT NAME | PACKAGING SIZE | ART. NO. |
|--|---------------------------|------------|
|  Cardboard Storage Box for 15 ml Centrifuge Tubes, 36 places Places: 36  Create PDF Data Sheet | 1 Piece (1 Box × 1 Piece) | E2023-1615 |
|  Cardboard Storage Box for 50 ml Centrifuge Tubes, 16 places Places: 16  Create PDF Data Sheet | 1 Piece (1 Box × 1 Piece) | E2023-1650 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk