




## Cardboard Storage Boxes for Centrifuge Tubes



### Product Highlights

- Moisture-resistant coating
- 16-place box for 50 ml tubes
- 36-place box for 15 ml tubes
- Removable dividers

## Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data

Art. No.	See variations
Dimensions (W × D × H)	146 mm × 146 mm × 76 mm
Dimensions (W × D)	146 mm × 146 mm
Max. tube length	140 mm
Autoclavable	No
Material	Coated Cardboard
Color	White
Sustainability attributes	Recyclable
For use with (SL)	Centrifuge Tubes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## More informations about Cardboard Storage Boxes for Centrifuge Tubes





White cardboard [storage boxes](#) with lift-off lid for the storage of 15 ml or 50 ml [centrifuge tubes](#). Moisture-resistant coating. Removable dividers. Box dimensions (WxDxH): 146 x 146 x 76 mm.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">Cardboard Storage Box for 15 ml Centrifuge Tubes, 36 places</a> <b>Places: 36</b>  <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Box × 1 Piece)	E2023-1615
 <a href="#">Cardboard Storage Box for 50 ml Centrifuge Tubes, 16 places</a> <b>Places: 16</b>  <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Box × 1 Piece)	E2023-1650

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)