




Cardboard Storage Boxes 133x133 Footprint, 50mm High



Product Highlights

- Moisture resistant coating
- Standard 133x133 footprint, for use in most freezer racking units
- Surface is easy to write on
- Removable dividers
- Recyclable

Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



General Data

Art. No.	See variations
Pack Size	1 Piece (1 Box × 1 Piece)
Dimensions (W × D × H)	133 mm × 133 mm × 50 mm
Dimensions (W × D)	133 mm × 133 mm
Autoclavable	No
Material	Coated Cardboard
Color	Various
Sustainability attributes	Made from responsibly sourced, biodegradable materials. Recyclable after use.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



More informations about Cardboard Storage Boxes 133x133 Footprint, 50mm High













Starlab's Coloured Cardboard Storage Boxes are made from high-quality responsibly sourced cardboard and coated with a moisture resistant coating, but are still easy to write on. The standard 133x133x50 footprint provides compatibility with most freezer racking. These boxes offer economical but safe and reliable storage at low temperatures. Available in various colours to aid colour-coding of your samples. Supplied pre-assembled, with removable grid insert. The boxes are fully recyclable after use.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk







All Variations

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	Cardboard Storage Box, 133x133 Footprint, 50mm High, 81-Place Grid, Yellow Places: 81 Color: Yellow  Create PDF Data Sheet	1 Piece (1 Box × 1 Piece)	E2023-8106
	Cardboard Storage Box, 133x133 Footprint, 50mm High, 81-Place Grid, Red Places: 81 Color: Red  Create PDF Data Sheet	1 Piece (1 Box × 1 Piece)	E2023-8104
	Cardboard Storage Box, 133x133 Footprint, 50mm High, 81-Place Grid, Green Places: 81 Color: Green  Create PDF Data Sheet	1 Piece (1 Box × 1 Piece)	E2023-8102
	Cardboard Storage Boxes 133x133 Footprint, 50mm High, 81-Place Grid, blue Places: 81 Color: Blue  Create PDF Data Sheet	1 Piece (1 Box × 1 Piece)	E2023-8101
	Cardboard Storage Box, 133x133 Footprint, 50mm High, 100-Place Grid, Yellow Places: 100 Color: Yellow  Create PDF Data Sheet	1 Piece (1 Box × 1 Piece)	E2023-3106
	Cardboard Storage Box, 133x133 Footprint, 50mm High, 100-Place Grid, Red Places: 100 Color: Red  Create PDF Data Sheet	1 Piece (1 Box × 1 Piece)	E2023-3104

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



PRODUCT NAME	PACKAGING SIZE	ART. NO.
 Cardboard Storage Box, 133x133 Footprint, 50mm High, 100-Place Grid, Green Places: 100 Color: Green  Create PDF Data Sheet	1 Piece (1 Box × 1 Piece)	E2023-3102
 Cardboard Storage Box, 133x133 Footprint, 50mm High, 100-Place Grid, Blue Places: 100 Color: Blue  Create PDF Data Sheet	1 Piece (1 Box × 1 Piece)	E2023-3101

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk