




## 10 ml Gilson-Style Graduated Macro Tip (non-sterile)



### Product Highlights

- From 1200 µl tips to 10 ml tips
- Made from polypropylene
- Compatible with a wide range of pipettes
- Different pack types available
- Certified

## Product Applications

Pipetting of large volumes of liquids

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## General Data

Art. No.	I1054-0001
Pack Size	250 Tips (10 Racks × 25 Tips)
Sterile	No
Autoclavable	Yes
Filter tip	Non-filter
Volume	10 ml
Material	Polypropylene
DNA free	Yes
DNase free	Yes
Pyrogen free	Yes
RNase free	Yes
Graduations	Yes
Design Tips	Gilson-Style
For use with	<b>Eppendorf</b> <sup>®</sup> : Reference 2, Research, Research plus, Xplorer <b>Gilson</b> <sup>®</sup> : Pipetman Concept, Pipetman P, Pipetman Ultra <b>HTL</b> : Discovery Comfort, Labmate Pro, Optipette <b>Rainin</b> <sup>®</sup> : Classic, Pipet-Lite, Pipet-Plus, Pipet-Plus with LTS <b>Socorex</b> <sup>®</sup> : Acura electro 955 <b>STARLAB</b> <sup>®</sup> : StarPet

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## More informations about 10 ml Gilson-Style Graduated Macro Tip (non-sterile)











Large Volume Tips

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)





## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">10 ml Finn/Socorex-Style Blue Tip (non-sterile)</a> <b>Design Tips:</b> Finn/Socorex-Style <b>Sterile:</b> No <b>Pack Size:</b> 100 Tips (1 Bag × 100 Tips)  <a href="#">Create PDF Data Sheet</a>	100 Tips (1 Bag × 100 Tips)	I1052-0000
 <a href="#">10 ml Finn/Socorex-Style Blue Tip (non-sterile)</a> <b>Design Tips:</b> Finn/Socorex-Style <b>Sterile:</b> No <b>Pack Size:</b> 250 Tips (10 Racks × 25 Tips)  <a href="#">Create PDF Data Sheet</a>	250 Tips (10 Racks × 25 Tips)	I1052-0001
 <a href="#">10 ml Finn/Socorex-Style Blue Tip (sterile)</a> <b>Design Tips:</b> Finn/Socorex-Style <b>Sterile:</b> Yes <b>Pack Size:</b> 250 Tips (10 Racks × 25 Tips)  <a href="#">Create PDF Data Sheet</a>	250 Tips (10 Racks × 25 Tips)	I1052-0010
 <a href="#">10 ml Gilson-Style Graduated Macro Tip (non-sterile)</a> <b>Design Tips:</b> Gilson-Style <b>Sterile:</b> No <b>Pack Size:</b> 100 Tips (1 Bag × 100 Tips)  <a href="#">Create PDF Data Sheet</a>	100 Tips (1 Bag × 100 Tips)	I1054-0000
 <a href="#">10 ml Gilson-Style Graduated Macro Tip (non-sterile)</a> <b>Design Tips:</b> Gilson-Style <b>Sterile:</b> No <b>Pack Size:</b> 250 Tips (10 Racks × 25 Tips)  <a href="#">Create PDF Data Sheet</a>	250 Tips (10 Racks × 25 Tips)	I1054-0001

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">10 ml Gilson-Style Graduated Macro Tip (sterile)</a> <b>Design Tips:</b> Gilson-Style <b>Sterile:</b> Yes <b>Pack Size:</b> 250 Tips (10 Racks × 25 Tips)  <a href="#">Create PDF Data Sheet</a>	250 Tips (10 Racks × 25 Tips)	I1054-0011

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)