




## 84 ml Reservoir, Multi-Well | 12-Channel | High Profile



### Product Highlights

- All robotic reservoirs are certified RNase, DNase and Pyrogen free

## Product Applications

PCR Setup

Assay-Setup

Serial Dilutions

Distributing Reagents

Sample or Reagent Transfer

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## General Data

Art. No.	E2999-8412
Pack Size	25 Reservoirs (1 Box × 25 Reservoirs)
Sterile	No
Volume	84 ml
Volume per trough	7 ml
Number of troughs	12
Material	Polypropylene
Bottoms	Trough-shape
Profile	High
Well	Multi-well

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)







## More informations about 84 ml Reservoir, Multi-Well | 12-Channel | High Profile

### High quality for high throughput!

Starlab reservoirs conform to the standard SBS footprint for compatibility with most robotic systems. They have virtually no dead volume and are available in high- or lowprofile and two well styles. The sturdy polypropylene construction ensures rigidity during automation and provides excellent chemical resistance.



## All Variations

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	<a href="#">240 ml Reservoir, Single-Well   96 Pyramid Bottoms   High Profile</a> <b>Bottoms:</b> Pyramid <b>Number of troughs:</b> 1   <a href="#">Create PDF Data Sheet</a>	25 Reservoirs (1 Box × 25 Reservoirs)	E2999-2496
	<a href="#">80 ml Reservoir, Multi-Well   8-Channel   High Profile</a> <b>Bottoms:</b> Trough-shape <b>Number of troughs:</b> 8   <a href="#">Create PDF Data Sheet</a>	25 Reservoirs (1 Box × 25 Reservoirs)	E2999-8008
	<a href="#">84 ml Reservoir, Multi-Well   12-Channel   High Profile</a> <b>Bottoms:</b> Trough-shape <b>Number of troughs:</b> 12   <a href="#">Create PDF Data Sheet</a>	25 Reservoirs (1 Box × 25 Reservoirs)	E2999-8412

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)