



## 100 – 1000 µl ErgoOne® Single-Channel Pipette



### Product Highlights

- Minimal operating forces
- Lightweight ergonomic design
- Colour coded
- Volume lock
- Universal tip fit (excl. 5.0 ml)
- Fully autoclavable
- 3 year warranty

## Product Applications

Pipetting liquids

Mixing liquids

Reverse pipetting of "tricky" liquids (e. g. with high viscosity)

Filling of plates, gels, and reaction vessels

Phase extraction and removal of supernatants

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## General Data

Art. No.	S7110-1000
Pipetting type	Air-cushion
Channel type	1-channel
Color code	Yellow
Operation	Mechanical, separate control button and ejector
Pack Size	1 Piece (1 Box × 1 Piece)
Autoclavable	Yes
Volume Setting	Variable
Cone type	Standard-cone
Volume range	100 - 1,000 µl
Volume	1000 µl

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## More informations about 100 – 1000 µl ErgoOne® Single-Channel Pipette

### Setting the new standard in pipetting

Every day, pipettes are confronted with multiple tasks that must be exact and precise. Reliability is key in the laboratory and your pipette must be trusted to perform well, day in, day out.

As one of the most frequently used instruments in the laboratory, pipettes must also be comfortable to use and provide the features you need. ErgoOne® pipettes combine researcher-requested features with the latest ergonomic principles to provide optimum performance.













ErgoOne® provides you with accuracy and precision you can trust, and a stylish new look you will love.

Lightweight and extremely easy to use, ErgoOne® can be trusted to provide accurate and precise pipetting all day, every day. Use with STARLAB's TipOne® Tips for the perfect pipetting partnership.

ErgoOne® is designed to fit all popular pipette tips (excl. 5.0 ml) without any loss of accuracy or precision.








## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">0.1 - 2.5µl ErgoOne® Single-Channel Pipette</a> <b>Volume range:</b> 0.1 - 2.5 µl <b>Cone type:</b> Standard-cone   <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Box × 1 Piece)	S7100-0125
 <a href="#">2 - 20µl ErgoOne® Single-Channel Pipette (Standard Cone)</a> <b>Volume range:</b> 2.0 - 20 µl <b>Cone type:</b> Standard-cone   <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Box × 1 Piece)	S7100-0220
 <a href="#">2 - 20µl ErgoOne® Single-Channel Pipette (Micro Cone)</a> <b>Volume range:</b> 2.0 - 20 µl <b>Cone type:</b> Micro-cone   <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Box × 1 Piece)	S7100-0221
 <a href="#">0.5 - 10µl ErgoOne® Single-Channel Pipette</a> <b>Volume range:</b> 0.5 - 10 µl <b>Cone type:</b> Standard-cone   <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Box × 1 Piece)	S7100-0510
 <a href="#">10 – 100 µl ErgoOne® Single-Channel Pipette</a> <b>Volume range:</b> 10 - 100 µl <b>Cone type:</b> Standard-cone   <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Box × 1 Piece)	S7100-1100
 <a href="#">20 – 200 µl ErgoOne® Single-Channel Pipette</a> <b>Volume range:</b> 20 - 200 µl <b>Cone type:</b> Standard-cone   <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Box × 1 Piece)	S7100-2200

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.


**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)








PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">100 – 1000 µl ErgoOne® Single-Channel Pipette</a> <b>Volume range:</b> 100 - 1,000 µl <b>Cone type:</b> Standard-cone  <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Box × 1 Piece)	S7110-1000
 <a href="#">500 – 5000 µl ErgoOne® Single-Channel Pipette</a> <b>Volume range:</b> 500 - 5,000 µl <b>Cone type:</b> Standard-cone  <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Box × 1 Piece)	S7150-5000

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## Accessories

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	<a href="#">ErgoOne® Wall Mount, Self-Adhesive, for one Pipette</a>	1 Pcs. (1 Pack × 1 Pcs.)	S7200-0001
	<a href="#">Self-Adhesive Shelf Holder for one ErgoOne®</a>	1 Piece (1 Pack × 1 Piece)	E3099-0001
	<a href="#">Benchtop Support for ErgoOne® Pipettes</a>	5 Pcs. (1 Pack × 5 Pcs.)	E3099-0002
	<a href="#">Acrylic Stand for up to 5 Pipettes (Universal)</a>	1 Piece (1 Pack × 1 Piece)	E3099-0000
	<a href="#">ErgoOne® Carousel Stand for up to 8 Pipettes (Universal)</a>	1 Piece (1 Box × 1 Piece)	S7200-0000

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)