



## 10 – 100 µl ErgoOne® 8-Channel Pipette



### Product Highlights

- Minimal operating forces
- Lightweight ergonomic design
- Colour coded
- Volume lock
- Universal tip fit
- Fully autoclavable
- 3 year warranty

## Product Applications

Pipetting liquids

Mixing liquids

Reverse pipetting of "tricky" liquids (e. g. with high viscosity)

Filling of plates, gels, and reaction vessels

Phase extraction and removal of supernatants

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## General Data

Art. No.	S7108-1100
Pack Size	1 Piece (1 Box × 1 Piece)
Color code	Yellow
Autoclavable	Yes
Volume range	10 - 100 µl
Volume	100 µl
Volume Setting	Variable
Channel type	8-channel
Operating mode	Manual
Pipetting type	Air-cushion
Operation	Mechanical, separate control button and ejector

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## More informations about 10 – 100 µl ErgoOne® 8-Channel Pipette

### A great tool to speed up your pipetting process

ErgoOne® multi-channel pipettes are designed in accordance with the latest ergonomic principles. The pipette is extremely easy to use and has minimal operating forces to help prevent hand fatigue. ErgoOne® multi-channel pipettes provide accurate and precise pipetting in a variety of volume ranges.













As one of the most frequently used instruments in the laboratory, pipettes must also be comfortable to use and provide the features you need. ErgoOne® pipettes combine researcher-requested features with the latest ergonomic principles to provide optimum performance. ErgoOne® gives you accuracy and precision you can trust, with a stylish new look you will love.

Lightweight and extremely easy to use, ErgoOne® can be trusted to provide accurate and precise pipetting all day, every day. Use with STARLAB's TipOne® pipette tips for the perfect pipetting partnership.

ErgoOne® is designed to fit all popular pipettes without any loss of accuracy and precision.



## All Variations





PRODUCT NAME		PACKAGING SIZE	ART. NO.
	<a href="#">0.5 – 10 µl ErgoOne® 8-Channel Pipette</a> <b>Volume range:</b> 0.5 - 10 µl <b>Channel type:</b> 8-channel   <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Box × 1 Piece)	S7108-0510
	<a href="#">10 – 100 µl ErgoOne® 8-Channel Pipette</a> <b>Volume range:</b> 10 - 100 µl <b>Channel type:</b> 8-channel   <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Box × 1 Piece)	S7108-1100
	<a href="#">30 – 300 µl ErgoOne® 8-Channel Pipette</a> <b>Volume range:</b> 30 - 300 µl <b>Channel type:</b> 8-channel   <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Box × 1 Piece)	S7108-3300
	<a href="#">0.5 – 10µl ErgoOne® 12-Channel Pipette</a> <b>Volume range:</b> 0.5 - 10 µl <b>Channel type:</b> 12-channel   <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Box × 1 Piece)	S7112-0510
	<a href="#">10 – 100µl ErgoOne® 12-Channel Pipette</a> <b>Volume range:</b> 10 - 100 µl <b>Channel type:</b> 12-channel   <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Box × 1 Piece)	S7112-1100
	<a href="#">30 – 300µl ErgoOne® 12-Channel Pipette</a> <b>Volume range:</b> 30 - 300 µl <b>Channel type:</b> 12-channel   <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Box × 1 Piece)	S7112-3300

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## Accessories

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	<a href="#">ErgoOne® Carousel Stand for up to 8 Pipettes (Universal)</a>	1 Piece (1 Box × 1 Piece)	S7200-0000
	<a href="#">Acrylic Stand for up to 5 Pipettes (Universal)</a>	1 Piece (1 Pack × 1 Piece)	E3099-0000
	<a href="#">Benchtop Support for ErgoOne® Pipettes</a>	5 Pcs. (1 Pack × 5 Pcs.)	E3099-0002
	<a href="#">Self-Adhesive Shelf Holder for one ErgoOne®</a>	1 Piece (1 Pack × 1 Piece)	E3099-0001

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)