



5.0 ml ErgoOne E/Finn-Style Pipette Tip (non-sterile), Loose



Product Highlights

- Specially designed for optimal tip fit on ErgoOne® E Pipettes
- Made from polypropylene
- Compatible with a wide range of pipettes
- Different pack types available
- Certified
- Compatible with ErgoOne E®, Brand® and Finnpipettes®

Product Applications

Pipetting of large volumes of liquids

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



General Data

Art. No.	E1070-2601
Sterile	No
Autoclavable	Yes
Filter tip	Non-filter
Volume	5 ml
Pack Size	1,000 Tips (5 Bags × 200 Tips)
Material	Polypropylene
Design Tips	Finn-Style
For use with brand + product (SL)	Abimed : Discovery Comfort, LABMATE Optima BRAND : Transferpette, Transferpette Electronic, Transferpette Pro, Transferpette S Eppendorf : Reference 2, Research pro, Xplorer Gilson : Pipetman Concept, Pipetman P, Pipetman Ultra HTL : Discovery Comfort, Labmate Pro, Optipette Rainin : Classic, EDP1, EDP-3 Plus, Pipet-Lite, Pipet-Plus STARLAB : ErgoOne E Thermo Fisher Scientific Inc. : Finn Focus, Finnpipette Digital, Finnpipette F1/F2, Finnpipette Novus

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



More informations about 5.0 ml ErgoOne E/Finn-Style Pipette Tip (non-sterile), Loose













Large Volume Tips

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk














All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <p>5 ml Biohit-Style Tip (non-sterile) Design Tips: Biohit-Style Sterile: No Pack Size: 250 Tips (1 Bag × 250 Tips)</p> <p> Create PDF Data Sheet</p>	250 Tips (1 Bag × 250 Tips)	E1009-7000
 <p>5 ml Biohit-Style-XL Tip (non-sterile) Design Tips: Biohit-Style-XL Sterile: No Pack Size: 250 Tips (1 Bag × 250 Tips)</p> <p> Create PDF Data Sheet</p>	250 Tips (1 Bag × 250 Tips)	E1009-7001
 <p>5 ml Finn-Style Tip (non-sterile) Design Tips: Finn-Style Sterile: No Pack Size: 250 Tips (1 Bag × 250 Tips)</p> <p> Create PDF Data Sheet</p>	250 Tips (1 Bag × 250 Tips)	I1009-5000
 <p>5 ml Gilson-Style Graduated Macro Tip (non-sterile) Design Tips: Gilson-Style Sterile: No Pack Size: 250 Tips (1 Bag × 250 Tips)</p> <p> Create PDF Data Sheet</p>	250 Tips (1 Bag × 250 Tips)	I1050-0000
 <p>5 ml Gilson-Style Graduated Macro Tip (non-sterile) Design Tips: Gilson-Style Sterile: No Pack Size: 500 Tips (10 Racks × 50 Tips)</p> <p> Create PDF Data Sheet</p>	500 Tips (10 Racks × 50 Tips)	I1050-0700
 <p>5 ml Gilson-Style Graduated Macro Tip (sterile) Design Tips: Gilson-Style Sterile: Yes Pack Size: 500 Tips (10 Racks × 50 Tips)</p> <p> Create PDF Data Sheet</p>	500 Tips (10 Racks × 50 Tips)	I1050-0710

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <p>5 ml Eppendorf-Style Graduated Macro Tip (non-sterile) Design Tips: Eppendorf-Style Sterile: No Pack Size: 250 Tips (1 Bag × 250 Tips)</p> <p> Create PDF Data Sheet</p>	250 Tips (1 Bag × 250 Tips)	I1053-0000
 <p>5 ml Eppendorf-Style Graduated Macro Tip (non-sterile) Design Tips: Eppendorf-Style Sterile: No Pack Size: 500 Tips (10 Racks × 50 Tips)</p> <p> Create PDF Data Sheet</p>	500 Tips (10 Racks × 50 Tips)	I1053-0700
 <p>5 ml Eppendorf-Style Graduated Macro Tip (sterile) Design Tips: Eppendorf-Style Sterile: Yes Pack Size: 500 Tips (10 Racks × 50 Tips)</p> <p> Create PDF Data Sheet</p>	500 Tips (10 Racks × 50 Tips)	I1053-0710
 <p>5.0 ml ErgoOne E/Finn-Style Pipette Tip (non-sterile), Loose Design Tips: Finn-Style Sterile: No Pack Size: 200 Tips (1 Bag × 200 Tips)</p> <p> Create PDF Data Sheet</p>	200 Tips (1 Bag × 200 Tips)	E1070-2596
 <p>5.0 ml ErgoOne E/Finn-Style Pipette Tip (non-sterile), Loose Design Tips: Finn-Style Sterile: No Pack Size: 1,000 Tips (5 Bags × 200 Tips)</p> <p> Create PDF Data Sheet</p>	1,000 Tips (5 Bags × 200 Tips)	E1070-2601
 <p>5 ml Gilson-Style Graduated Macro Tip (non-sterile) Design Tips: Gilson-Style Sterile: No Pack Size: 2,500 Tips (10 Bag × 250 Tips)</p> <p> Create PDF Data Sheet</p>	2,500 Tips (10 Bag × 250 Tips)	I1050-0000-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk