




5 ml Gilson-Style Graduated Macro Filter Tip (sterile)



Product Highlights

- From 1200 µl tips to 10 ml tips
- Made from polypropylene
- Compatible with a wide range of pipettes
- Different pack types available
- Certified
- DNase free, RNase free, DNA free, Pyrogen free

Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



General Data

Art. No.	I1050-0718
Pack Size	500 Tips (10 Racks × 50 Tips)
Sterile	Yes
Filter tip	Filter
Volume	5 ml
Material	Polypropylene
DNA free	Yes
DNase free	Yes
Human DNA free	Yes
Pyrogen free	Yes
RNase free	Yes
Graduations	No
Design Tips	Gilson-Style
For use with	Abimed: Discovery Comfort, LABMATE Optima Eppendorf: Reference 2, Research plus, Research pro, Xplorer Gilson: Pipetman Concept, Pipetman P, Pipetman Ultra HTL: Discovery Comfort, Labmate Pro, Optipette Nichiryo: 31N, NPLO, NPX Rainin: Classic, EDP1, EDP-3 Plus, Pipet-Lite, Pipet-Plus Socorex: Acura electro 935 STARLAB: ErgoOne, StarPet Thermo Fisher Scientific Inc.: Finn Focus, Finnpipette F1/F2

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



More informations about 5 ml Gilson-Style Graduated Macro Filter Tip (sterile)











Large Volume Filter Tips

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk





All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 5 ml Gilson-Style Graduated Macro Filter Tip (sterile) Design Tips: Gilson-Style Sterile: Yes Pack Size: 500 Tips (10 Racks × 50 Tips)  Create PDF Data Sheet	500 Tips (10 Racks × 50 Tips)	I1050-0718
 5 ml Eppendorf-Style Graduated Macro Filter Tip (sterile) Design Tips: Eppendorf-Style Sterile: Yes Pack Size: 500 Tips (10 Racks × 50 Tips)  Create PDF Data Sheet	500 Tips (10 Racks × 50 Tips)	I1053-0718
 5 ml Eppendorf-Style Graduated Macro Filter Tip (sterile) Design Tips: Eppendorf-Style Sterile: Yes Pack Size: 100 Tips (2 Racks × 50 Tips)  Create PDF Data Sheet	100 Tips (2 Racks × 50 Tips)	I1053-0728
 5 ml Gilson-Style Graduated Macro Filter Tip (sterile) Design Tips: Gilson-Style Sterile: Yes Pack Size: 100 Tips (2 Racks × 50 Tips)  Create PDF Data Sheet	100 Tips (2 Racks × 50 Tips)	I1050-0728
 5.0 ml ErgoOne E/Finn-Style Filter Tip (sterile), Racked Design Tips: Finn-Style Sterile: Yes Pack Size: 50 Tips (1 Rack × 50 Tips)  Create PDF Data Sheet	50 Tips (1 Rack × 50 Tips)	E1009-5718

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



PRODUCT NAME	PACKAGING SIZE	ART. NO.
 10 ml Gilson-Style Graduated Macro Filter Tip (sterile) Design Tips: Gilson-Style Sterile: Yes Pack Size: 50 Tips (2 Racks × 25 Tips)  Create PDF Data Sheet	50 Tips (2 Racks × 25 Tips)	I1054-0028

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk