



ErgoOne® Starter Pack 4 (10)



Product Highlights

- Four ErgoOne® single-channel pipettes: **10µl, 20µl (standard), 200µl & 1000µl**
- One rack of 10/20ul XL Graduated TipOne® Tips, (Sterile)
- One rack of 200µl UltraPoint Graduated TipOne® Tips, (Sterile)
- One rack of 1250 µl XL Graduated TipOne® Tips, (Sterile)
- Four Self-Adhesive Acrylic Pipette Shelf Holders
- One STARLAB Lab Timer

Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



General Data

Art. No.	STARTERPACK410
Volume range	0.5 - 1,000 µl
Volume	1000 µl
Volume Setting	Variable
Pack Size	1 Piece (1 Box × 1 Piece)
Channel type	1-channel
Operating mode	Manual

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



More informations about ErgoOne® Starter Pack 4 (10)

Every day, pipettes are confronted with multiple tasks that must be exact and precise. Reliability is key in the laboratory and your pipette must be trusted to perform well, day in, day out. ErgoOne® pipettes provide:

- › Minimal operating forces
- › Lightweight ergonomic design
- › Colour coded
- › Volume lock
- › Universal tip fit







Start your new pipetting experience with ErgoOne® single-channel pipettes and original TipOne® tips - the perfect pipetting partnership. All supplied in an eco-friendly shipper box made with 30% grass! This material requires less water and energy consumption and significantly less chemical use.

ErgoOne® Starter Pack 4 (10) contains:

- 1x S7100-0510 (0.5 - 10µl ErgoOne® Single Channel Pipette)
- 1x S7100-0220 (2.0 - 20µl ErgoOne® Single Channel Pipette, standard cone)
- 1x S7100-2200 (20 - 200µl ErgoOne® Single Channel Pipette)
- 1x S7110-1000 (100 - 1000µl ErgoOne® Single Channel Pipette)
- 3x racks of sterile TipOne tips
- 4x E3099-0001 (Self-Adhesive Shelf Holder for 1 Pipette)
- 1x STARLAB Lab Timer



All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 ErgoOne® Starter Pack 4 (2.5) Volume range: 0.1 - 1,000 µl  Create PDF Data Sheet	1 Piece (1 Box × 1 Piece)	STARTERPACK42.5
 ErgoOne® Starter Pack 4 (10) Volume range: 0.5 - 1,000 µl  Create PDF Data Sheet	1 Piece (1 Box × 1 Piece)	STARTERPACK410
 ErgoOne® Starter Pack Volume range: 20 - 1,000 µl  Create PDF Data Sheet	1 Piece (1 Box × 1 Piece)	STARTERPACK

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk