



## 5 ml Eppendorf-Style Graduated Macro Filter Tip (sterile)



### Product Highlights

- From 1200 µl tips to 10 ml tips
- Made from polypropylene
- Compatible with a wide range of pipettes
- Different pack types available
- Certified
- free of DNase , RNase , DNA and Pyrogene
- DNase Testing is by gel electrophoresis. Test sensitivity 10<sup>-7</sup> Kunitz units/µl.
- RNase Testing is by gel electrophoresis. Test sensitivity 10<sup>-9</sup> Kunitz units/µl.
- DNA Testing is by gel electrophoresis. Test sensitivity 30 pg.
- Pyrogens Testing is for endotoxins by the LAL (Limulus Amebocyte Lysate)
- Gel Cloth method. Test sensitivity 0.06 EU/ml (Endotoxin Units).

## Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## General Data

Art. No.	I1053-0728
Pack Size	100 Tips (2 Racks × 50 Tips)
Sterile	Yes
Filter tip	Filter
Volume	5 ml
Material	Polypropylene
DNA free	Yes
DNase free	Yes
Human DNA free	Yes
Pyrogen free	Yes
RNase free	Yes
Graduations	Yes
Design Tips	Eppendorf-Style
For use with	<b>Abimed</b> : Discovery Comfort, LABMATE Optima <b>Biohit</b> : eLine, ePet, mLine, Picus, Proline <b>Eppendorf</b> : Reference 2, Research plus, Research pro, Xplorer <b>Gilson</b> : Pipetman Concept, Pipetman P, Pipetman Ultra <b>Hamilton</b> : SofTouch <b>HTL</b> : Discovery Comfort, Labmate Pro, Optipette <b>Nichiryo</b> : 31N, NPLO, NPX <b>Rainin</b> : Classic, EDP1, EDP-3 Plus, Pipet-Lite, Pipet-Plus <b>Sartorius</b> : eLine, ePet, mLine, Picus, Picus NxT, Proline, Proline Plus, Tacta <b>Socorex</b> : Acura 835F

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## More informations about 5 ml Eppendorf-Style Graduated Macro Filter Tip (sterile)











Large Volume Filter Tips

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)





## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">5 ml Gilson-Style Graduated Macro Filter Tip (sterile)</a> <b>Design Tips:</b> Gilson-Style <b>Sterile:</b> Yes <b>Pack Size:</b> 500 Tips (10 Racks × 50 Tips)  <a href="#">Create PDF Data Sheet</a>	500 Tips (10 Racks × 50 Tips)	I1050-0718
 <a href="#">5 ml Eppendorf-Style Graduated Macro Filter Tip (sterile)</a> <b>Design Tips:</b> Eppendorf-Style <b>Sterile:</b> Yes <b>Pack Size:</b> 500 Tips (10 Racks × 50 Tips)  <a href="#">Create PDF Data Sheet</a>	500 Tips (10 Racks × 50 Tips)	I1053-0718
 <a href="#">5 ml Eppendorf-Style Graduated Macro Filter Tip (sterile)</a> <b>Design Tips:</b> Eppendorf-Style <b>Sterile:</b> Yes <b>Pack Size:</b> 100 Tips (2 Racks × 50 Tips)  <a href="#">Create PDF Data Sheet</a>	100 Tips (2 Racks × 50 Tips)	I1053-0728
 <a href="#">5 ml Gilson-Style Graduated Macro Filter Tip (sterile)</a> <b>Design Tips:</b> Gilson-Style <b>Sterile:</b> Yes <b>Pack Size:</b> 100 Tips (2 Racks × 50 Tips)  <a href="#">Create PDF Data Sheet</a>	100 Tips (2 Racks × 50 Tips)	I1050-0728
 <a href="#">5.0 ml ErgoOne E/Finn-Style Filter Tip (sterile), Racked</a> <b>Design Tips:</b> Finn-Style <b>Sterile:</b> Yes <b>Pack Size:</b> 50 Tips (1 Rack × 50 Tips)  <a href="#">Create PDF Data Sheet</a>	50 Tips (1 Rack × 50 Tips)	E1009-5718

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">10 ml Gilson-Style Graduated Macro Filter Tip (sterile)</a> <b>Design Tips:</b> Gilson-Style <b>Sterile:</b> Yes <b>Pack Size:</b> 50 Tips (2 Racks × 25 Tips)  <a href="#">Create PDF Data Sheet</a>	50 Tips (2 Racks × 25 Tips)	I1054-0028

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)