



## ErgoOne® FAST Benchtop Support



### Product Highlights

## Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data

Art. No.	P7166-6306
Pack Size	1 Piece (1 Bag x 1 Piece)
For use with (SL)	ErgoOne® FAST Pipette Controller

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## More informations about ErgoOne® FAST Benchtop Support

During development of the [ErgoOne](#)® FAST electronic pipette controller, a high pipetting speed was recognised as an important feature. The result is convincing! The [ErgoOne](#)® FAST can fill a 25 ml pipette in less than 4 seconds!

Every aspect of [ErgoOne](#)® FAST has been developed with ergonomics in mind. [ErgoOne](#)® FAST is one of the lightest pipette controllers on the market. It is designed for comfort and balances perfectly in the hand to provide fatigue-free pipetting, even during long periods.

This benchtop support fits perfectly with the [ErgoOne](#)® FAST.