



## 2 ml Serological Pipette (sterile), individually wrapped



### Product Highlights

- Accuracy  $\pm 2\%$
- Bidirectional graduations (excl. 2.0 ml) provide added applicability
- Non-cytotoxic and non-haemolytic USP Class VI polystyrene
- Gamma sterilized to SAL  $10^{-6}$  (ISO11137)
- Certified RNase, DNase, latex and pyrogen free
- Colour-coded for easy identification
- Negative graduations (excl. 100 ml) allows additional working volume
- Supplied with a polyolefin filter plug
- Individually wrapped in traceable paper-plastic peel-back packs showing lot number and expiry date
- Specially designed carton with perforated tear-strip doubles as efficient, bench top dispenser
- ISO 9001 and ISO 13485 manufactured

## Product Applications

Serial dispensing of aliquots of different volumes

Resuspension of bacteria or cell pellets and nucleic acid precipitates

Aspiration of cell layer from, e.g., Ficoll® gradient

Phase extraction and removal of lower phases

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data

Art. No.	E4860-0002
Length	272.0 mm
Sterile	Yes
Volume range	1 - 2 ml
Volume	2 ml
Pack Size	500 Pcs. (1 Pack × 500 Pcs.)
Color code	Green
Material	Polystyrene
DNase free	Yes
Pyrogen free	Yes
RNase free	Yes
Graduations	Yes
Color	Clear
For use with (SL)	ErgoOne® FAST Pipette Controller

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## More informations about 2 ml Serological Pipette (sterile), individually wrapped

### Discover the possibilities!















STARLAB's new extended range of Serological Pipettes cover seven different volumes and, for the first time, includes a 100 ml pipette. The pipettes are constructed from polystyrene for sterile serological and tissue culture applications. Colour-coding on both the pipettes and the packaging provides easy identification. The black printed, easy to read ascending and descending graduations simplify fluid measuring and dispensing, and all pipettes (excl. 100 ml) feature negative graduations for extra working volume. Perfect for use with all popular pipette controllers, including ErgoOne® FAST.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">1 ml Serological Pipette (sterile), individually wrapped</a> <b>Volume:</b> 1 ml  <a href="#">Create PDF Data Sheet</a>	500 Pcs. (1 Pack × 500 Pcs.)	E4860-0001
 <a href="#">2 ml Serological Pipette (sterile), individually wrapped</a> <b>Volume:</b> 2 ml  <a href="#">Create PDF Data Sheet</a>	500 Pcs. (1 Pack × 500 Pcs.)	E4860-0002
 <a href="#">5 ml Serological Pipette (sterile), individually wrapped</a> <b>Volume:</b> 5 ml  <a href="#">Create PDF Data Sheet</a>	200 Pcs. (1 Pack × 200 Pcs.)	E4860-0005
 <a href="#">10 ml Serological Pipette (sterile), individually wrapped</a> <b>Volume:</b> 10 ml  <a href="#">Create PDF Data Sheet</a>	200 Pcs. (1 Pack × 200 Pcs.)	E4860-0010
 <a href="#">25 ml Serological Pipette (sterile), individually wrapped</a> <b>Volume:</b> 25 ml  <a href="#">Create PDF Data Sheet</a>	150 Pcs. (1 Pack × 150 Pcs.)	E4860-0025
 <a href="#">50 ml Serological Pipette (sterile), individually wrapped</a> <b>Volume:</b> 50 ml  <a href="#">Create PDF Data Sheet</a>	100 Pcs. (1 Pack × 100 Pcs.)	E4860-0050
 <a href="#">100 ml Serological Pipette (sterile), individually wrapped</a> <b>Volume:</b> 100 ml  <a href="#">Create PDF Data Sheet</a>	50 Pcs. (1 Pack × 50 Pcs.)	E4860-1100

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)