



5.0 ml Aspiration Pipettes (sterile), Ind. wrapped



Product Highlights

- Non-cytotoxic and non-haemolytic USP Class VI polystyrene
- Gamma sterilized to SAL 10^{-6} (ISO 11137)
- Certified RNase, DNase, latex and pyrogen free
- Individually wrapped in traceable paper-plastic peel-back packs showing lot number and expiry date
- Specially designed carton with perforated tear-strip doubles as efficient, bench top dispenser
- ISO 9001 and ISO 13485 manufactured

Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



General Data

| | |
|-------------------|---|
| Art. No. | E4861-0005 |
| Sterile | Yes |
| Volume range | 2 - 5.0 ml |
| Volume | 5.0 ml |
| Pack Size | 200 Individually wrapped sachets (1 Box × 200 Individually wrapped sachets) |
| Material | Polystyrene |
| For use with (SL) | Aspirator |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



More informations about 5.0 ml Aspiration Pipettes (sterile), Ind. wrapped







Sterile, polystyrene pipette without a filter or graduations. Aspiration pipettes are used for vacuum aspirating procedures or transfers of fluids where accuracy is not a requirement.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



All Variations

| PRODUCT NAME | PACKAGING SIZE | ART. NO. |
|--|--|------------|
|  2.0 ml Aspiration Pipette (sterile), Individually wrapped Volume: 2.0 ml  Create PDF Data Sheet | 500 Individually wrapped sachets (1 Box × 500 Individually wrapped sachets) | E4861-0002 |
|  5.0 ml Aspiration Pipettes (sterile), Individually wrapped Volume: 5.0 ml  Create PDF Data Sheet | 200 Individually wrapped sachets (1 Box × 200 Individually wrapped sachets) | E4861-0005 |
|  10.0 ml Aspiration Pipettes (sterile), Individually wrapped Volume: 10.0 ml  Create PDF Data Sheet | 200 Individually wrapped sachets (1 Box × 200 Individually wrapped sachets) | E4861-0010 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk