



## Sealing ring for rotor lid of rotor "Standard 24 AT"



### Product Highlights

- Spare part for use with Rotor Lid of the Microcentrifuge 24 Rotor Standard 24 AT
- Fully autoclavable

## Product Applications

Sample

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## General Data

|           |                           |
|-----------|---------------------------|
| Art. No.  | S8030-0012                |
| Pack Size | 1 Piece (1 Box × 1 Piece) |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## More informations about Sealing ring for rotor lid of rotor "Standard 24 AT"

Microcentrifuge 24 Rotor Standard 24 AT replacement part.

Content: 2x sealing rings, 1 x lubricant tube, and instructions.

An illustrated description of how to replace the sealing ring on the rotor lid can also be found in the download section.