




Microplate Mixer



Product Highlights

- Suitable for use with all skirted plates with SBS footprint
- Wide range of speed from 0 to 1500 rpm
- Small footprint
- Stepless speed control
- Various accessories for microtitre plates and tubes
- 2-year warranty

Product Applications

Sample preparation

Sample mixing

Resuspending cell pellets

Resuspending nucleic acid pellets

Resuspending chemical pellets

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



General Data

Art. No.	N2400-8040
Dimensions (W × D × H)	150 mm × 260 mm × 80 mm
Dimensions (W × D)	150 mm × 260 mm
Orbital Diameter	4.5 mm
Weight	3 kg
Pack Size	1 Piece (1 Box × 1 Piece)
Component Type	Orbital
Power supply	100 – 240 V, 50/60 Hz
Max. power consumption	20 W
Operating temperature	4 to 40 °C
Permissible Ambient Temperature	4 – 40 °C
Max. Load	0.5 kg
Max. RCF	1500 rpm
Speed Range	0 – 1,500 rpm

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



More informations about Microplate Mixer

Compact and versatile!






This microplate mixer is ideal for the mixing and shaking of any skirted plate (PCR, [Cell Culture](#), Deepwell, Microtitre and [Microplates](#)). The mixer can hold either one or two plates, depending on the attachment being used, and can also mix tubes with the special tube adapter accessories. The stepless, variable vibration speed is excellent for applications such as the resuspending of pellets. The small footprint saves space on your bench.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



Accessories

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 Attachment for two microtiterplates	1 Piece (1 Box × 1 Piece)	N2400-8042
 Microplate Mixer Adapter Plate	1 Piece (1 Box × 1 Piece)	N2400-8043
 Tube Adapter for 15 x 1.5/ 2.0ml Tubes	1 Piece (1 Box × 1 Piece)	N2400-6012
 Tube Adapter for 8 x 15ml Centrifuge Tubes	1 Piece (1 Box × 1 Piece)	N2400-6013
 Attachment for one microtiterplate	1 Piece (1 Box × 1 Piece)	N2400-8041

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk