



## Mini Dry Bath Metal Block for 6 x 15ml Tubes



### Product Highlights

## Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data

|                        |                            |
|------------------------|----------------------------|
| Art. No.               | N2400-4025                 |
| Dimensions (W × D × H) | 71 mm × 47 mm × 75 mm      |
| Dimensions (W × D)     | 71 mm × 47 mm              |
| Well diameter          | 17.3 mm                    |
| Weight                 | 0.6 kg                     |
| Well depth             | 70 mm                      |
| Pack Size              | 1 Piece (1 Box × 1 Piece)  |
| System                 | Mini Dry Bath Incubator    |
| For use with (SL)      | 6 x 15 ml centrifuge tubes |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## More informations about Mini Dry Bath Metal Block for 6 x 15ml Tubes

### Small, but powerful!

A small but powerful unit suitable for a range of applications. The digital microprocessor controls the preset temperature of the metal block and the fast ramping of temperature saves your valuable time. Supplied with a transparent lid for temperature uniformity, the unit can also be used as a water bath.

We offer a range of accessories below.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)