



## Hand controller for Aspirator



### Product Highlights

- Ergonomic design
- Fully autoclavable
- Flow through suction tube - media is not in contact with the interior of the hand controller

## Product Applications

Aspiration of Liquids

Cell-Culture

Sterile Filtration

De-gassing of buffers

Solidphase Extraction

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## General Data

Art. No.	N2400-9001
Pack Size	1 Piece (1 Box × 1 Piece)

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## More informations about Hand controller for Aspirator


Hand controller for Aspirator N2400-9000. To connect a second hand controller to the Aspirator, use either connection kit N2400-9006 or Quick coupling N2400-9008. Includes Adapter for pipettes, adapter for pipette tips, wall mount for hand controller and 2.5 m tubing.

The adapter can be used:

- For pipettes (Pasteur, aspiration, and graduated; Outer Diameter 5.5–8.0 mm)
- For pipette tips (2–200 µL and 50–1000 µL)



## Accessories

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">Adapter for pipettes</a>	1 Piece (1 Box × 1 Piece)	N2400-9016

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)