




200 µl Round Gel-Loading Tip, Racked (non-sterile)



Product Highlights

- Flexible, long and narrow
- Loading samples onto gels (eg. polyacrylamide or agarose)
- Handling of small volumes
- Compatible with a wide range of pipettes
- Different pack types available
- Certified

Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



General Data

Art. No.	I1022-0600
Length	68 mm
Sterile	No
Filter tip	Non-filter
Volume	200 µl
Pack Size	1,020 Tips (5 Racks × 204 Tips)
Tip Type	round gel-loading
DNA free	Yes
DNase free	Yes
Pyrogen free	Yes
RNase free	Yes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



More informations about 200 µl Round Gel-Loading Tip, Racked (non-sterile)









All tips are certified RNase and DNase free. Packaging of these tips differs from the [TipOne](#)[®] range.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 200 µl Round Gel-Loading Tip, Loose (non-sterile) Sterile: No Pack Size: 1,000 Tips (1 Bags × 1000 Tips)  Create PDF Data Sheet	1,000 Tips (1 Bags × 1000 Tips)	I1022-0000
 200 µl Round Gel-Loading Tip, Loose (non-sterile), Case Sterile: No Pack Size: 10,000 Tips (10 Bags × 1000 Tips)  Create PDF Data Sheet	10,000 Tips (10 Bags × 1000 Tips)	I1022-0000-C
 200 µl Round Gel-Loading Tip, Racked (non-sterile) Sterile: No Pack Size: 1,020 Tips (5 Racks × 204 Tips)  Create PDF Data Sheet	1,020 Tips (5 Racks × 204 Tips)	I1022-0600
 200 µl Round Gel-Loading Tip, Racked (Sterile) Sterile: Yes Pack Size: 1,020 Tips (5 Racks × 204 Tips)  Create PDF Data Sheet	1,020 Tips (5 Racks × 204 Tips)	I1022-0610

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk