



## Underbench Freezer-86°C



### Product Highlights

- Use as your personal -80°C storage place
- Eco-friendly, HCFC-free Nature R refrigerant
- USB port for temperature logger
- No filter; maintenance free!
- 48 hour battery backup
- 2-year warranty
- High / low temperature alarm
- Porthole for external probe & CO2 backup
- Spark-free motor
- Supplied with 2 adjustable shelves
- 91 liters Volume

## Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data

Art. No.	N2700-1000
Dimensions (W × D × H)	595 mm × 690 mm × 825 mm
Dimensions (W × D)	595 mm × 690 mm
Weight	60 kg
Pack Size	1 Piece (1 Box × 1 Piece)
Power supply	220 V
Max. temp.	-86 °C
Temperature range	-60 to -86 °C
Temperature min. (SL)	-60 °C
Output	220V–240V, 50/60 Hz

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## More informations about Underbench Freezer-86°C

Keep it chilled!

At last, a small -86°C freezer that is affordable!

Store your samples in this space saving, eco-friendly unit directly under your lab bench. The freezer is compatible with a range of storage boxes and racks. Limited maintenance, whilst easy to install and move. A great choice if you are looking for a small size, energy efficient freezer with a great storage capacity.

The freezer is supplied with two adjustable shelves, or use one of our stainless steel racks that fit perfectly inside the freezer to utilise storage capacity to the full. A rack is available to hold 100 mm high storage boxes with standard 133 x 133 mm footprint.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## Accessories

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">U100 Rack for 100 mm High Boxes, 33 places</a>	1 Piece (1 Box × 1 Piece)	N2700-1003

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)