



## Calibration on Purchased Timer



### Product Highlights

- Calibration of Timer (C0003-0000) must be ordered at the same time as Timer (N9900-0001), on the same Purchase Order No.
- Your new Timer (N9900-0000) will then be calibrated prior to shipping
- Certificate of Calibration included

## Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## General Data

Art. No.	C0003-0000
Pack Size	1 Piece (1 Box × 1 Piece)
Dimensions (W × D × H)	73 mm × 13 mm × 59 mm
Dimensions (W × D)	73 mm × 13 mm
Operating modes	4
Timer	Count down / Count up
Clock	12 hr/24 hr

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## More informations about Calibration on Purchased Timer

### Timing is everything!

Timer/Clock with a large easy to read LED display. There are four channels on which cycles can be independently or simultaneously set. You can easily toggle between your pre-set channels and the clock in order to keep on top of your experiments, or to know when it's lunch!

### Need a calibrated timer?

To receive a pre-calibrated Timer: Calibration of Timer (C0003-0000) must be ordered at the same time as Timer (N9900-0001), on the same Purchase Order No. We will then know to calibrate the Timer prior to shipping, and send it with a Certificate of Calibration.