



## StarStore 49



### Product Highlights

- 49-place polypropylene storage box with fixed 7x7 dividers for tubes up to 16mm
- Transparent lift-off lid with frosted writing panel on side
- Vent holes in base to aid draining and lid removal
- Optional acetate grid matrix available to sit under the lid
- Temperature use: -90 °C to 121 °C
- Autoclavable
- Dimensions: 130 x 130 x 45mm
- Max. tube height: 60 mm

## Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## General Data

Art. No.	E2349-5009
Pack Size	5 Pcs. (1 Pack × 5 Pcs.)
Dimensions (W × D × H)	130 mm × 130 mm × 45 mm
Dimensions (W × D)	130 mm × 130 mm
Max. tube length	60 mm
External width	130 mm
Autoclavable	Yes
Material	Polypropylene
Color	Natural
Places	49

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## More informations about StarStore 49

49-place polypropylene storage box with fixed, 7x7 dividers for tubes up to 16mm. The non-hinged, lift-off lid enables tubes of varying heights to be stored inside. Standard 130x130mm footprint fits standard freezer racking.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)