



StarStore 100-Place Storage Box, Natural



Product Highlights

- 100-place polypropylene storage box with fixed 10x10 dividers for tubes up to 13mm
- Hinged, lockable lid with positioning lugs for stable stacking
- Moulded reference grid on base.
- Temperature use -90°C to 121°C
- Autoclavable
- Dimensions (WxDxH): 140 x 155 x 55 mm

Product Applications

Sample storage

Storage of microcentrifuge tubes or screw cap microtubes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



General Data

Art. No.	I2310-5800
Dimensions (W × D × H)	140 mm × 155 mm × 55 mm
Dimensions (W × D)	140 mm × 155 mm
Autoclavable	Yes
Pack Size	5 Pcs. (1 Pack × 5 Pcs.)
Material	Polypropylene
Color	Natural
Places	100

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



More informations about StarStore 100-Place Storage Box, Natural















100-place polypropylene storage box with fixed, 10x10 dividers for tubes up to 13mm, eg. 1.5 / 2.0 ml microcentrifuge tubes or screw cap microtubes. Hinged, lockable lid with positioning lugs for stable stacking. Can be used in the freezer racking made specifically for use with this storage box.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk



All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 StarStore 100-Place Storage Box, Natural Color: Natural  Create PDF Data Sheet	5 Pcs. (1 Pack × 5 Pcs.)	I2310-5800
 StarStore 100-Place Storage Box, Neon Blue Color: Neon blue  Create PDF Data Sheet	5 Pcs. (1 Pack × 5 Pcs.)	I2310-5841
 StarStore 100-Place Storage Box, Neon Green Color: Neon green  Create PDF Data Sheet	5 Pcs. (1 Pack × 5 Pcs.)	I2310-5842
 StarStore 100-Place Storage Box, Neon Orange Color: Neon orange  Create PDF Data Sheet	5 Pcs. (1 Pack × 5 Pcs.)	I2310-5843
 StarStore 100-Place Storage Box, Neon Pink Color: Neon pink  Create PDF Data Sheet	5 Pcs. (1 Pack × 5 Pcs.)	I2310-5844
 StarStore 100-Place Storage Box, Neon Yellow Color: Neon yellow  Create PDF Data Sheet	5 Pcs. (1 Pack × 5 Pcs.)	I2310-5846
 StarStore 100-Place Storage Box, Neon Mixed* Color: Neon mixed  Create PDF Data Sheet	5 Pcs. (1 Pack × 5 Pcs.)	I2310-5848

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: infoline@starlab.co.uk